

TL-232P

Trial Lens Set



MAGNON

TL-232P Trial Lens Set

MAGNON

The MAGNON TL-232P is double curve trial lens set which offers a standard selection of lens powers to meet practitioner's need.

FEATURES

- Double curve lenses of 22mm in diameter
- Color coded plastic rims of 38mm in diameter
- Removable plastic inner tray with space for trial frame
- Aluminum carrying case with combination locks

SPECIFICATIONS		
Spheres	34 pairs each of plus and minus	
Cylinders	19 pairs each of plus and minus	
Prisms	10 pieces	
Accessories	10 pieces	
	● : Occluder2	G : Green lens 1
	○ : Pin hole (1mm)1	F : Frosted lens 1
	○ : Pin hole (1.5mm)1	N : 0.0D lens 1
	⊓ : Slit (1.5mm)1	M : White Maddox 1
	R : Red lens1	

- Spheres	-Cylinders	+ Cylinders	+ Spheres
20	6.00	6.00	20
18	5.00	5.00	18
16	4.50	4.50	16
15	4.00	4.00	15
14	3.50	3.50	14
13	3.25	3.25	13
12	3.00	3.00	12
11	2.75	2.75	11
10	2.50	2.50	10
9.00	2.25	2.25	9.00
8.00	2.00	2.00	8.00
7.00	1.75	1.75	7.00
6.50	1.50	1.50	6.50
6.00	1.25	1.25	6.00
5.50	1.00	1.00	5.50
5.00	0.75	0.75	5.00
4.50	0.50	0.50	4.50
4.00	0.25	0.25	4.00
3.75	0.12	0.12	3.75
3.50			3.50
3.25			3.25
3.00	Accessories	Prisms	3.00
2.75	● ● ○ ○ ⊓ R G F N M	1° 2° 3° 4° 5° 6° 7° 8° 9° 10°	2.75
2.50			2.50
2.25			2.25
2.00			2.00
1.75			1.75
1.50			1.50
1.25			1.25
1.00			1.00
0.75			0.75
0.50			0.50
0.25			0.25
0.12			0.12

Specifications and design are subject to change for improvement without notice.

H. OGINO & CO., LTD.

Yokohama Port P. O.Box 170, Yokohama 231-8691, Japan
 Tel: (045)662-5261 Fax: (045)651-5467
 E-mail : magnon@h-ogino.co.jp URL : www.h-ogino.co.jp