





TECNIS™ IOL with OptiEdge™ Design
Hydrophobic Acrylic Posterior Chamber
Model ZA9003, 6 mm Optic, Mod C Haptic

Model ZA9003





DESCRIPTION		ZA9003
OPTIC CHARACTERISTICS		
Powers:		+10.0 to +30.0 Diopters in .5 increments
Diameter:		6 mm
Shape:		Biconvex, anterior modified prolate
Material:		UV blocking hydrophobic acrylic
Manipulation Holes:		None
A-Constant ¹ :		119.1
Theoretical AC Depth ¹ :		5.6 mm
Refractive Index:		1.47
Surgeon Factor ² :		1.85
HAPTIC CHARACTERISTICS		
Overall Length:		13 mm
Style:		Modified C
Material:		60% Blue core Polymethylmethacrylate (PMMA) Monofilament
Angle:		5°
Compressibility (10 mm) ³ :		173 mg
Average Weight in Air ¹ :		21.5 mg
RECOMMENDED INSERTION INSTRUMENTS		MODEL
THE UNFOLDER™ Emerald Series Handpiece		EMERALDT & EMERALDXL
THE UNFOLDER™ Emerald Series Cartridge		EMERALDC



To order or for more information,
call your AMO representative or
visit www.TECNISIOL.com.

1. Value theoretically derived for a typical 20.0 D lens.
2. Calculated based on Holladay I formula (Holladay JT, Prager TC, Chandler TY, Musgrove KH, Lewis JW, Ruiz RS. A three-part system for referring intraocular lens power calculations. *J Cataract Refract Surg.* 1988;14(1):17-24).
3. Data on file, Advanced Medical Optics, Inc.

The TECNIS™ intraocular lens (ZA9003) is produced under at least one of the following U.S. Letters Patent Nos.: 5,104,590; 5,444,106; 6,007,747; 6,609,793; and 6,705,729. Patents pending. Foreign equivalents available upon request.

The AMO ADVANCED MEDICAL OPTICS INNOVATING VISION logo, THE UNFOLDER and TECNIS are trademarks of Advanced Medical Optics, Inc. OptiEdge is a trademark of Ocular Sciences, Inc. Rx Only. ©2005 ADVANCED MEDICAL OPTICS, INC., Santa Ana, CA 92705 www.amo-inc.com

Johnson & Johnson

AMO
ADVANCED MEDICAL OPTICS
Innovating Vision