

FORCEPS



CILIA.....	101
CONJUNCTIVA.....	102
SUTURING.....	103
CORNEAL.....	104
TYING.....	105
FIXATION.....	109
IRIS.....	110
CAPSULORHEXIS.....	113
IOL.....	121
MUSCLE.....	126
VITREORETINAL.....	129
CHALAZION.....	134
MOSQUITO.....	135
SERREFINES · TOWEL CLAMPS.....	136

G-18720 *HENRY*
CILIA FORCEPS
straight jaws
smooth, blunt, squared tips
overall length 100 mm



G-18730 *GRADLE*
CILIA FORCEPS
curved jaws
smooth, blunt, squared tips
overall length 98 mm



G-18750 *BARRAQUER*
CILIA FORCEPS
lightly angled
6.5 mm platform
overall length 109 mm



G-18740 *BEER*
CILIA FORCEPS
3 mm platform
rounded tips 4.5 mm wide
overall length 78 mm



G-18850 *SHAAF*
CILIA AND FOREIGN BODY FORCEPS
grooved tips
curved and pointed
finger grip handle
overall length 96 mm



G-31620 *BARRAQUER*
CILIA AND SUTURE FORCEPS
grooved tips, 7.0 mm platform
curved and pointed
finger grip handle
overall length 112 mm




G-32740 *HAMBURG MODEL*
CILIA AND FOREIGN BODY FORCEPS
tapered fine tips
overall length 90 mm




G-31230

KREFELD MODEL
CONJUNCTIVA FORCEPS
 standard
 1 x 2 teeth
 overall length 100 mm




G-31235

FREIBURG MODEL
CONJUNCTIVA FORCEPS
 delicate
 1 x 2 teeth
 overall length 90 mm




G-18847

KAUFMANN
CONJUNCTIVA FORCEPS
 standard
 1 x 2 teeth
 5 mm tying platform
 overall length 75mm




G-18810

DRAEGER
CONJUNCTIVA FORCEPS
 1,2 mm wide
 2 atraumatic fixation knobs
 overall length 97 mm



G-18848

FOERSTER
CONJUNCTIVA FORCEPS
 delicate
 1 x 2 teeth
 finger grip handle
 overall length 95 mm






CASTROVIEJO
SUTURING FORCEPS

1 x 2 teeth
6 mm tying platform
overall length 100 mm


- G-19001** 0.28 mm, teeth 0.12 mm
- G-19005** 0.3 mm
- G-19010** 0.6 mm
- G-19015** 0.9 mm
- G-19020** 1.2 mm




CASTROVIEJO
SUTURING FORCEPS

new model
1 x 2 teeth
6 mm tying platform
overall length 109 mm

- G-19022** 0.25 mm, teeth 0.12 mm
- G-19023** 0.3 mm
- G-19024** 0.6 mm
- G-19025** 0.9 mm
- G-19026** 1.2 mm



- G-19008** *MUNICH MODEL*
SUTURING FORCEPS
1 x 2 teeth, 0.35 mm
5.5 mm tapered tying platform
overall length 100 mm



ZUERICH MODEL
SUTURING FORCEPS

delicate
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 75 mm
straight
G-18965
G-18966 angled



G-18950 COLIBRI FORCEPS
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 75 mm



G-18960 MICRO COLIBRI FORCEPS
delicate
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 77 mm



10

SUTURING · COLIBRI STYLE · CORNEAL FIXATION

G-15134 MICRO COLIBRI FORCEPS
extra delicate
1 x 2 teeth, 0.12 mm
4 mm tying platform
overall length 76 mm



G-18975 BARRAQUER SUTURE REMOVING FORCEPS
very fine pointed tips
4 mm tying platform
overall length 77 mm



G-18972 THOMAS NEUHANN KERATOPLASTY FORCEPS
for corneal suturing and fixation
two prongs with 1 x 2 teeth
0.8 mm prong spread
overall length 80mm



G-18970

REMKY

KERATOPLASTY FORCEPS

two prongs with 1 x 1 teeth, angled
1.7 mm prong spread
for corneal suturing and fixation
overall length 80 mm




G-38623

REMKY

DMEK KERATOPLASTY FORCEPS

modified by *SZURMAN*
two prongs with 1 x 1 teeth, straight
1.7 mm prong spread
33 mm shafts
to hold graft during preparation
overall length 89 mm



G-18940


PAUFIQUE

UNIVERSAL SUTURING FORCEPS

1 x 2 teeth
6 mm tying platform
overall length 88 mm

G-18945

standard
delicate




G-19045

PAUFIQUE


TYING FORCEPS

with platform
overall length 88 mm
standard, 6 mm platform
G-19046 delicate, 5 mm platform






G-19030 *TUEBINGEN MODEL*
TYING FORCEPS
 straight, smooth jaws
 6 mm tying platform
 overall length 101 mm



G-19040 *TUEBINGEN MODEL*
TYING FORCEPS
 curved, smooth jaws
 6 mm tying platform
 overall length 99 mm




10


TYING




G-19035 *TUEBINGEN MODEL*
TYING FORCEPS
 delicate, smooth jaws
 4 mm tying platform
 overall length 84mm




G-19034 *BONN MODEL*
TYING FORCEPS
 extra delicate, smooth jaws
 3 mm tying platform
 overall length 73 mm



G-19032 **TYING FORCEPS**
angled, 7.5 mm
6 mm tying platform
overall length 99 mm




G-31625 *DARDENNE*
TYING FORCEPS
angled, 6.5 mm
6 mm tying platform
round handle with guide pin
overall length 102 mm



G-19038 *BANGERTER*
TYING FORCEPS
9.5 mm platform
longitudinal v-shape groove
overall length 100 mm



G-19048 **TYING FORCEPS**
utility serrated jaws
1 mm wide
overall length 106 mm



G-32947

HATTENBACH
HYBRID TYING FORCEPS
 angled, for the anterior chamber
 intraocular 1.5 mm micro tying tips
 25 gauge / 0.5 mm
 overall length 98 mm



for use in combination with
 hybrid needle holder G-32949

G-31623

KELMAN-MCPHERSON
FRAGMENT FORCEPS
 angled, smooth jaws
 7.5 mm platform
 overall length 100 mm



G-33941

KELMAN-MCPHERSON
FRAGMENT FORCEPS
 angled, fine grooved jaws,
 for femto-LASIK to retrieve lenticule
 during SMILE procedures
 7.5 mm platform
 overall length 100 mm



G-31622

GASKIN
FRAGMENT FORCEPS
 angled
 with tying platform
 9 mm platform, overall length 107 mm
 G-31624 11 mm platform, overall length 100 mm



G-37700 **JEWELERS TYPE FORCEPS**
standard, straight
overall length 120 mm



G-37705 **JEWELERS TYPE FORCEPS**
delicate, straight
overall length 115 mm



G-37710 **JEWELERS TYPE FORCEPS**
extra delicate, straight
overall length 115 mm



G-19080 *BARRAQUER*
FIXATION FORCEPS
3 sets of 1 x 2 teeth
11 mm fixation platform
overall length 100 mm



GRAEF
FIXATION FORCEPS
4,5 mm wide
8 interlocking teeth
overall length 110 mm
G-19050 without lock
G-19055 with lock




G-19060 *HEIDELBERG MODEL*
FIXATION FORCEPS
3,5 mm wide
6 interlocking teeth
overall length 110 mm



DRESSING FORCEPS

serrated jaws
overall length 100 mm


- G-18780 0.5 mm
- G-18781 0.7 mm
- G-18782 1.0 mm
- G-18783 1.2 mm
- G-18784 1.4 mm



TISSUE FORCEPS

1 x 2 teeth
overall length 100 mm


- G-18790 0.5 mm
- G-18791 0.7 mm
- G-18792 1.0 mm
- G-18793 1.2 mm
- G-18794 1.4 mm



IRIS FORCEPS

serrated, delicate
overall length 73 mm
straight
curved

- G-18855
- G-18856




BONN MODEL

IRIS FORCEPS

serrated, extra delicate
overall length 73 mm
straight
angled at 6 mm
angled at 8 mm

- G-18908
- G-18906
- G-18911




BONN MODEL

IRIS FORCEPS

1 x 2 teeth, 0.12 mm
extra delicate
overall length 73mm
straight
angled at 6 mm
angled at 8 mm

- G-18900
- G-18905
- G-18910




BONN MODEL

IRIS FORCEPS

long, straight
1 x 2 teeth, 0.12 mm
overall length 102 mm
extra delicate
delicate, with 4.5 mm tying platform

G-18845

G-18890




G-18840

SWISS MODEL

IRIS FORCEPS

curved, 1 x 2 teeth
overall length 100 mm




G-18860

HEIDELBERG MODEL

IRIS FORCEPS

light curve, 1 x 2 teeth
overall length 72 mm




G-18870

HEIDELBERG MODEL

IRIS FORCEPS


strong curve, 1 x 2 teeth
overall length 70 mm



G-18857

IRIS FORCEPS


straight, delicate
1 x 2 teeth
overall length 72 mm




G-31953

IRIS FORCEPS


straight, strong
1 x 2 teeth
overall length 72 mm




G-18846 *HAMBURG MODEL*
IRIS FORCEPS
 straight, delicate
 1 x 2 teeth
 overall length 103 mm



G-32240 *SWISS MODEL*
IRIS FORCEPS
 straight, extra delicate
 1 x 2 teeth
 overall length 70 mm




G-18880 *HESS-ENGELKING*
IRIS FORCEPS
 1 x 2 teeth
 overall length 72 mm




10

IRIS FIXATION

G-18897 *ELZE*
IRIS FORCEPS
 extra delicate
 1 x 2 teeth
 overall length 78 mm




G-18913 *THOMAS NEUHANN*
IRIS FORCEPS
 curved shaft
 1 x 2 teeth
 23 gauge / 0.6 mm
 overall length 95 mm



G-31307

CAPSULORHEXIS FORCEPS


strong, with heel stop to protect the iris
interlocking micro serrated tips
12 mm straight shafts
overall length 103 mm



G-31308

CAPSULORHEXIS FORCEPS

round handle with guide pin
extra delicate with
interlocking micro serrated tips
13 mm curved shafts
overall length 105 mm




G-31300

MASKET

CAPSULORHEXIS FORCEPS

delicate, with heel stop to protect the iris
interlocking micro serrated tips
12 mm curved shafts
overall length 103 mm



G-31304


CAPSULORHEXIS FORCEPS

new model, with guide pin
extra delicate with
interlocking micro serrated tips
13 mm curved shafts
overall length 103 mm







G-31301 CAPSULORHEXIS FORCEPS
 extra delicate
 with smooth triangular grasping tips
 13 mm curved shafts
 overall length 98 mm




G-31305 CAPSULORHEXIS FORCEPS
 extra delicate with
 interlocking micro serrated tips
 13 mm curved shafts
 overall length 98 mm



G-31303 H.-R. KOCH CAPSULORHEXIS FORCEPS
 extra delicate with
 interlocking micro serrated tips
 13 mm curved shafts
 1 mm marks from tip to bent
 overall length 98 mm



G-31299 CAPSULORHEXIS FORCEPS
 short handle, extra delicate
 cross serrated grasping tips
 12 mm straight shafts
 overall length 83 mm





G-32926

AUFFARTH

MICS DOUBLE CROSS ACTION CAPSULORHEXIS FORCEPS

interlocking micro serrated tips
12 mm curved shafts, angled 130°
with overlap system
prevents inadvertent iris capture
for 1.6 mm incisions
overall length 100 mm




G-32928

INAMURA

MICS CAPSULORHEXIS FORCEPS

cross action
interlocking micro serrated tips
10.5 mm curved shafts tip to pivot
angled 130°
overall length 103 mm




G-31306

KAMPIK

MICS CAPSULORHEXIS FORCEPS

22 gauge / 0.7 mm curved shaft
90° cystotome tips, reverse action
overall length 93 mm




G-32935

H.-R. KOCH


MICS CAPSULE FORCEPS

short handle
22 gauge / 0.7 mm curved shaft with
micro grasping tips
overall length 98 mm







G-32938 *H.-R. KOCH*
MICS CAPSULE FORCEPS
 22 gauge / 0.7 mm curved shaft with
 micro grasping tips
 6 mm titanium handle, blue-colored
 overall length 144 mm



G-32936 *TOBIAS NEUHANN*
MICRO FORCEPS
 multifunction asymmetrical tips with
 grasping platform
 23 gauge / 0.6 mm straight shaft
 6 mm titanium handle, blue-colored
 overall length 155 mm



G-32937 *IKEDA*
MICS CAPSULORHEXIS FORCEPS
 23 gauge / 0.6 mm curved shaft with
 micro grasping tips
 6 mm titanium handle, blue-colored
 overall length 144 mm




IKEDA
CAPSULORHEXIS FORCEPS
 interchangeable tips
 23 gauge / 0.6 mm curved shaft
 overall length 40 mm

G-38244 textured, micrograsping tips "BEAK"
G-38242 textured, micrograsping tips



G-38244
 "beak"



G-38242



for use with handle G-38240

G-38240 *JAPAN MODEL*
HANDLE
 for capsule scissors and forceps
 interchangeable tips
 6 mm, titanium
 overall length 90 mm




G-32944

HATTENBACH

CAPSULORHEXIS FORCEPS

25 gauge / 0.5 mm angled shaft with
micro grasping tips
overall length 88 mm




G-18805

TOBIAS NEUHANN

DESCEMET MEMBRANE FORCEPS

delicate blunt grasping tips
with platform
overall length 75 mm




G-38624

SZURMAN

DMEK PREPARATION FORCEPS

for removal of Descemet layer
with smooth round tips to prevent tearing
the donor endothelium
overall length 75 mm




G-38628

YÖRÜK

DMEK TRANSPLANT FORCEPS

atraumatic grasping and removal of
Descemet's membrane
grasping area 6 mm
overall length 77 mm



G-38634

GERTEN


DSAEK/DMEK RHEXIS FORCEPS

for Descemetorhexis
reversed triangular tips to
score and remove the DM from the host cornea
23 gauge / 0.6 mm curved shaft
overall length 142 mm






G-38603 **DSAEK IMPLANTATION FORCEPS**
 13 mm curved shafts with limited closing
 to avoid injury of endothelial transplant
 1.5 mm serrated platform
 overall length 107 mm



G-38606 *THOMAS NEUHANN*
DSAEK FORCEPS
 1 x 2 teeth, 0.12 mm
 21 gauge / 0.8 mm curved shaft
 for endothelial transplant procedures
 overall length 97 mm




10


DSAEK · LASIK FLAP



G-33428 *PFAEFFL*
FLAPORHEXIS FORCEPS
 2.5 mm x 1.5 mm cross serrated tips, oval
 to lift femto-LASIK flaps
 overall length 103 mm




G-33429 *KNORZ*
FEMTO-LASIK FLAP FORCEPS
 curved 3.25 mm T-bar
 textured
 to grasp and lift LASIK flaps
 overall length 102 mm






G-31899 **LASIK CORNEAL FLAP FORCEPS**
angled shafts, 9 mm
highly polished round and smooth shafts
overall length 108 mm




G-18914 *BURATTO*
LASIK CORNEAL FLAP FORCEPS
curved shafts
cross serrated tips, 0.8 mm diameter
overall length 100 mm



G-33960 *RUPAL SHAH*
LENTICULE FORCEPS
3.5 mm serrated tips
20 gauge / 0.9 mm curved shaft
to retrieve lenticule during SMILE procedures
6 mm titanium handle, blue-colored
overall length 144 mm




G-33961 **LENTICULE FORCEPS**
7.5 mm shafts
with fine grooves
angled to side
to retrieve lenticule during SMILE procedures
overall length 99 mm







G-31518 **NUCLEUS SPREADING FORCEPS**
 cross-action
 3.0 mm textured paddle tips
 for cracking the nucleus
 overall length 112 mm




G-31519 *AKAHOSHI*
NUCLEUS SPREADING FORCEPS
 phaco prechopper
 4.5 mm delicate sharp tips, 1 mm profile
 2.0 mm cross-action
 for cracking, dividing, and rotating the nucleus
 overall length 105 mm



G-31720 *AKAHOSHI*
PRECHOPPER
 for micro-incisions, straight shafts, cross action
 round handle
 overall length 125 mm



G-31722 *AKAHOSHI*
HYBRID COMBO PRECHOPPER
 for micro-incisions, cross action, straight
 round handle
 length of tip 5.3 mm
 overall length 120 mm







G-31724 *AKAHOSHI*
UNIVERSAL PRECHOPPER
for micro-incisions, cross action, straight
round handle
overall length 120 mm




G-33440 *VON BARSEWISCH*
ANTERIOR CHAMBER FORCEPS
extra delicate
8.5 mm angled shafts
1.5 mm platform
overall length 72 mm



G-32165 *FECHNER*
IMPLANTATION FORCEPS
smooth curved shafts
overall length 80 mm



G-18935 *DODEN*
IMPLANTATION FORCEPS
curved
delicate lens holding shafts
overall length 110 mm







G-18930 *SHEPARD*
IMPLANTATION FORCEPS
 curved, 1.4 mm disc-shaped tips
 delicate lens holding shafts
 overall length 120 mm




G-31962 *T. KOHNEN*
IMPLANTATION FORCEPS
 silicone tube coated
 10 mm angled shafts
 overall length 107 mm



G-18918 *DICK*
FIXATION AND HOLDING FORCEPS
 for Artisan/Verisyse IOL
 with guide mark for lens maximum hold
 and stability
 18 gauge / 1.2 mm curved shaft
 overall length 143



G-18919 *DICK*
IOL HAPTIC EXPANDER
 for enclavation of iris-fixated IOLs
 21 gauge / 0.8 mm curved shaft
 overall length 140 mm





10


IOL IMPLANTATION · SILICONE COATED




G-18944 *T. NEUHANN*
IMPLANTATION FORCEPS
for retropupillary implantation of Iris-Claw-IOL
0.9 mm wide holding surface
20 gauge / 0.9 mm curved shaft
overall length 143 mm



G-18928 *SEKUNDO*
IMPLANTATION FORCEPS
for retropupillary implantation of Iris-Claw-IOL
20 gauge / 0.9 mm curved shaft with
lens holding plate, three-point fixation
and IOL centration guide mark
6 mm titanium handle, blue-colored
overall length 143 mm



G-31902 **FOLDER**
foldable IOL insertion forceps
cross action
6.5 mm concave jaws
overall length 103 mm




G-31903 *SPALECK*
FOLDER
foldable IOL insertion forceps
for silicone and acrylic IOLs
angled shafts, concave jaws
overall length 102 mm





G-31908 **FOLDER**
IOL folding forceps
for easy pickup and folding acrylic IOLs
overall length 105 mm



G-31994 *H.-R. KOCH*
PREFOLDER
for foldable silicone, hydrogel or acrylic IOLs
overall length 135 mm



10


IOL FOLDERS · IOL HOLDERS



G-31905 *H.-R. KOCH*
FOLDER
for foldable 3-piece IOLs
6,0 mm concave jaws
overall length 92 mm




G-31901 **HOLDING FORCEPS**
for foldable IOLs
highly polished round and smooth shafts
overall length 108 mm






G-18941 *STROBEL*
HOLDING FORCEPS
for Pharmacia CeeOn Edge IOL
overall length 110 mm




G-32823 **EXPLANATION FORCEPS**
for soft IOLs
2.0 mm x 1.4 mm paddle tip with
interlocking serrations and guide pin
overall length 108 mm



G-32997 *SNYDER-OSHER*
IOL FORCEPS
to hold the IOL while cutting it inside
the capsular bag
2 mm long jaws with 0.5 mm serrated tips
22 gauge / 0.7 mm lightly curved shaft
overall length 99 mm




G-18820 *REMKY*
SUPERIOR RECTUS FORCEPS
angled
1 x 2 teeth
overall length 109 mm




G-18800

ELSCHNIG
SUPERIOR RECTUS FORCEPS
 standard, straight
 1 x 2 teeth
 overall length 110 mm



G-18795

ELSCHNIG
SUPERIOR RECTUS FORCEPS
 delicate, straight
 1 x 2 teeth
 overall length 100 mm




10

SUPERIOR RECTUS · MUSCLE


G-19105

LEYDHECKER
SUPERIOR RECTUS FORCEPS
 curved, tip width 0.8 mm
 1 x 2 teeth, 1.1 mm
 overall length 110 mm



G-32335

CASTROVIEJO
SUPERIOR RECTUS FORCEPS
 angled with tying platform
 1 x 2 teeth, 0.9 mm
 overall length 100 mm



G-19112

BANGERTER
MUSCLE FORCEPS
 10 teeth
 self retaining
 overall length 100 mm



SATTLER

MUSCLE FORCEPS

with slide lock
overall length 100 mm
G-19120 3 teeth, bulbous tip, right
G-19125 3 teeth, bulbous tip, left



G-32350

TROELTSCH

NASAL FORCEPS

serrated tips
for packing and dressing during ENT procedures
overall length 150 mm



G-19107

ESSEN MODEL

CERCLAGE FORCEPS

strongly curved
3 mm platform with dual grooves
overall length 104 mm



G-33752

WATZKE

SPREADING FORCEPS

to stretch silicone sleeves
cross action handle and serrated tips
overall length 115 mm




G-33434

PFENNIGWERTH

HOLDING FORCEPS


for vitreoretinal lenses
curved with concave lens holding tips
overall length 110 mm



G-33439

SCLERAL PLUG FORCEPS


cross action handle, angled jaws to securely hold plugs during insertion and removal
overall length 105 mm



G-31570

SCLERAL PLUG FORCEPS

to securely hold plugs during insertion and removal
overall length 100 mm



10


RETINAL (SCLERAL PLUG, RETINAL TACK)



G-33432

SCLERAL PLUG


serrated,
20 gauge / 0.9 mm



G-33436

SCLERAL PLUG

to suture onto sclera
20 gauge / 0.9 mm



G-33435

SCLERAL PLUG

20 gauge / 0.9 mm



G-33441

SCLERAL PLUG

23 gauge / 0.6 mm

G-33438

FORCEPS

for HEIMANN retinal tacks
overall length 100 mm



G-33437T

**HEIMANN
RETINAL TACK**


titanium
- for single use only -




G-37420 *BARTZ-SCHMIDT*
TROCAR CLOSURE FORCEPS
with function for bulbus fixation
overall length 72 mm



G-37430 **TROCAR FORCEPS**
for trocar systems with groove
overall length 110 mm




G-33831 *REIMANN*
VITREORETINAL FORCEPS
for foreign body removal
with retractable magnet and 3-prongs
14 gauge / 2.1 mm straight shaft
shaft length 35 mm
overall length 138 mm




20 GAUGE

VITREORETINAL FORCEPS


20 gauge / 0.9 mm
 6 mm and 8 mm diameter handles
 titanium, blue-colored
 overall length 155 mm



straight
 serrated tips
G-36202 6 mm handle




straight
 textured tips
G-36206 6 mm handle
G-36208 8 mm handle




straight
 serrated "crocodile" tips
G-36256 6 mm handle
G-36258 8 mm handle

10


VITREORETINAL 20 GAUGE




straight
 silicone coated
G-36212 6 mm handle
G-36214 8 mm handle




TANO
 straight, serrated tips
G-36356 6 mm handle
G-36358 8 mm handle




A. MUELLER
 for retinal vessels / sheatotomy
G-36332 6 mm handle





straight
 grip tips
G-36216 6 mm handle
G-36218 8 mm handle





90° angle
 pointed tips
G-36276 6 mm handle
G-36278 8 mm handle





NAITO
 ILM
G-36316 6 mm handle
G-36318 8 mm handle





straight
 endgripping, extra delicate
 0.2 mm x 0.05 mm
G-36292 6 mm handle
G-36294 8 mm handle



straight
 endgripping, extra delicate
 0.3 mm x 0.05 mm
G-36296 6 mm handle
G-36298 8 mm handle



straight
 endgripping, delicate
 0.2 mm x 0.1 mm
G-36302 6 mm handle
G-36304 8 mm handle




straight
 endgripping
 0.3 mm x 0.1 mm
G-36306 6 mm handle
G-36308 8 mm handle


20 GAUGE

VITREORETINAL FORCEPS


20 gauge / 0.9 mm
6 mm and 8 mm diameter handles
titanium, blue-colored
overall length 155 mm




90° angle
membrane peeler
G-36226 6 mm handle



45° angle
membrane peeler
G-36232 6 mm handle




45° angle
G-36222 6 mm handle




membrane peeler tips
G-36252 6 mm handle
G-36254 8 mm handle




curved
membrane dissector
G-36266 6 mm handle



45° angle
membrane dissector
G-36272 6 mm handle




45° angle
horizontal subretinal
G-36236 6 mm handle




45° angle, vertical tips
subretinal
G-36246 6 mm handle
G-36248 8 mm handle




angled, extra delicate
3.0 mm tips
G-36336 6 mm handle



angled, extra delicate
4.5 mm tips
G-36326 6 mm handle
G-36328 8 mm handle



straight
circular tips
G-36242 6 mm handle




foreign body forceps
G-36262 6 mm handle
G-36264 8 mm handle

23 GAUGE


VITREORETINAL FORCEPS

23 gauge / 0.6 mm
6 mm and 8 mm diameter handles
titanium, blue-colored
overall length 155 mm




straight
serrated "crocodile" tips

G-36342 6 mm handle
G-36344 8 mm handle





TANO
straight, serrated tips

G-36362 6 mm handle
G-36364 8 mm handle




straight
endgripping, extra delicate
0.2 mm x 0.05 mm

G-36282 6 mm handle
G-36284 8 mm handle



straight
grip tips


G-36286 6 mm handle
G-36288 8 mm handle



25 GAUGE


VITREORETINAL FORCEPS

25 gauge / 0.5 mm
6 mm and 8 mm diameter handles
titanium, blue-colored
overall length 155 mm




straight
serrated "crocodile" tips

G-36322 6 mm handle
G-36324 8 mm handle




TANO
straight, serrated tips

G-36366 6 mm handle
G-36368 8 mm handle



straight
endgripping, extra delicate
0.2 mm x 0.05 mm

G-36312 6 mm handle
G-36314 8 mm handle




ERGOLINE

ERGOLINE FORCEPS


360° rotatable and interchangeable tips
overall length 55 mm
for Ergoline handpieces 6 mm and 8 mm diameters

ERGOLINE HANDPIECE



for Ergoline forceps
blue-colored
G-36600 6 mm titanium handpiece
G-36601 8 mm titanium handpiece





ERGOLINE FORCEPS
straight
serrated "crocodile" tips
G-36655 20 gauge / 0.9 mm
G-36740 23 gauge / 0.6 mm
G-36720 25 gauge / 0.5 mm





TANO
ERGOLINE FORCEPS
straight, serrated tips
G-36755 20 gauge / 0.9 mm
G-36765 23 gauge / 0.6 mm





NAITO
ERGOLINE FORCEPS ILM
G-36715 20 gauge / 0.9 mm





ERGOLINE FORCEPS
endgripping
extra delicate 0.2 mm x 0.05 mm
G-36690 20 gauge / 0.9 mm
G-36680 23 gauge / 0.6 mm
G-36710 25 gauge / 0.5 mm



ERGOLINE FORCEPS
endgripping
extra delicate 0.3 mm x 0.05 mm
G-36695 20 gauge / 0.9 mm



ERGOLINE FORCEPS
endgripping
delicate 0.2 mm x 0.1 mm
G-36700 20 gauge / 0.9 mm



ERGOLINE FORCEPS
endgripping
delicate 0.3 mm x 0.1 mm
G-36705 20 gauge / 0.9 mm

**SINGLE-USE
INTERCHANGEABLE FORCEPS
AVAILABLE ON P. 219**



G-18680 *LAMBERT*
CHALAZION FORCEPS
 round, 9 mm inside diameter
 solid and fenestrated plates
 overall length 92 mm



G-18670 *AYER*
CHALAZION FORCEPS
 round, 8 mm inside diameter
 interlacing tips
 solid and fenestrated plates
 overall length 90 mm



DESMARRES
CHALAZION FORCEPS
 oval
 solid and fenestrated plates
G-18685 20 mm, overall length 88 mm
G-18690 25 mm, overall length 90 mm
G-18695 31 mm, overall length 95 mm



10

CHALAZION • PTOSIS

G-18700 *CAUER*
CHALAZION FORCEPS
 oval, 13 mm x 11 mm inside diameter
 interlacing serrated tips
 solid and fenestrated plates
 overall length 99 mm



G-18710 *LORDAN*
CHALAZION FORCEPS
 triangular, 13 mm x 8 mm inside diameter
 solid and fenestrated plates
 overall length 95 mm



G-18715 *BERKE*
PTOSIS FORCEPS
 22 mm long serrated tips
 slide lock mechanism
 overall length 110 mm



HARTMANN
MOSQUITO FORCEPS

serrated
18 mm jaws
overall length 100 mm
straight
curved

G-18170
G-18171



HALSTEAD
MOSQUITO FORCEPS

serrated
25 mm jaws
overall length 125 mm
straight
curved

G-18180
G-18181



HALSTEAD
MOSQUITO FORCEPS

1 x 2 teeth
serrated
25 mm jaws
overall length 125 mm
straight
curved

G-18185
G-18182





G-18160 *BACKHAUS*
TOWEL FORCEPS
overall length 90 mm



G-33200 *SCHAEDEL*
TOWEL FORCEPS
overall length 90 mm

10

TOWEL CLAMPS · SERREFINES



G-18190 *DIEFFENBACH*
SERREFINE
straight, overall length 35 mm



G-18195 *DIEFFENBACH*
SERREFINE
curved, overall length 35 mm



G-31420 *DIEFFENBACH*
SERREFINE
straight, overall length 50 mm



G-31415 *DIEFFENBACH*
SERREFINE
curved, overall length 50 mm



G-18197 **SERREFINE**
extra delicate, overall length 30 mm