

NEEDLE HOLDERS



BARRAQUER

NEEDLE HOLDER

10 mm delicate jaws
0.8 x 0.55 mm tip
overall length 120 mm
curved, without lock
curved, with lock

G-17500
G-17501



BARRAQUER

NEEDLE HOLDER

10 mm delicate jaws
0.8 x 0.45 mm tip
overall length 120 mm
curved, without lock
curved, with lock

G-17505
G-17506



JACOBI

NEEDLE HOLDER

5 mm extra delicate jaws
1.0 x 0.5 mm tips
overall length 115 mm
straight, without lock
straight, with lock

G-17451
G-17453



JACOBI

NEEDLE HOLDER

5 mm extra delicate jaws
1.0 x 0.5 mm tips
overall length 115 mm
curved, without lock
curved, with lock

G-17450
G-17452



BARRAQUER

NEEDLE HOLDER

7.5 mm extra delicate jaws
0.6 x 0.3 mm tip
overall length 105 mm
straight, without lock
straight, with lock

G-17461
G-17463



BARRAQUER

NEEDLE HOLDER

7.5 mm extra delicate jaws
0.6 x 0.3 mm tip
overall length 105 mm
curved, without lock
curved, with lock

G-17460
G-17462



BARRAQUER-TROUTMAN
NEEDLE HOLDER

10 mm extra delicate jaws
0.8 x 0.6 mm tip
handle diameter 6 mm
overall length 100 mm
straight, without lock
straight, with lock

G-32260
G-32261



BARRAQUER-TROUTMAN
NEEDLE HOLDER

10 mm extra delicate jaws
0.8 x 0.6 mm tip
handle diameter 6 mm
overall length 100 mm
curved, without lock
curved, with lock

G-17535
G-17530



BARRAQUER
NEEDLE HOLDER

11 mm delicate jaws
0.9 x 0.8 mm tips
handle diameter 8 mm
overall length 140 mm
curved, without lock
curved, with lock

G-17510
G-17511



BARRAQUER
NEEDLE HOLDER

11 mm extra delicate jaws
0.9 x 0.55 mm tips
handle diameter 8 mm
overall length 140 mm
curved, without lock
curved, with lock

G-17515
G-17516



BARRAQUER
NEEDLE HOLDER
11 mm delicate jaws
1.1 x 0.75 mm tips
handle diameter 10 mm
overall length 140 mm
curved, without lock
curved, with lock

G-17520
G-17525



BARRAQUER
NEEDLE HOLDER
curved, without lock
8 mm standard jaws
1.2 x 0.8 mm tips
original model with longitudinal groove
handle diameter 10 mm
overall length 130 mm



G-32315



BARRAQUER
NEEDLE HOLDER
12 mm standard jaws
1.6 x 1.0 mm tips
handle diameter 10 mm
overall length 140 mm
curved, without lock
curved, with lock

G-32320
G-32321

CASTROVIEJO
NEEDLE HOLDER

short model
9 mm extra delicate jaws
0.9 x 0.55 mm tips
overall length 90 mm
straight, without lock
straight, with lock

G-17545
G-17540



CASTROVIEJO
NEEDLE HOLDER

10.5 mm delicate jaws
1.5 x 0.8 mm tips
overall length 140 mm
straight, without lock
straight, with lock

G-17560
G-17565




CASTROVIEJO
NEEDLE HOLDER

10.5 mm delicate jaws
1.5 x 0.8 mm tips
overall length 140 mm
curved, without lock
curved, with lock

G-17570
G-17575




COMBINED NEEDLE HOLDER / SCISSORS
curved 12 mm blades with
3 mm needle holder and tying platform
used to hold, place and cut fine suture material
overall length 125 mm
without lock
with lock



G-17440
G-17441



HATTENBACH
HYBRID NEEDLE HOLDER
angled, for the anterior chamber
25 gauge / 0.5 mm
overall length 90 mm



G-32949

for use in combination with
hybrid tying forceps G-32947



KALT
NEEDLE HOLDER
standard jaws
serrated, with longitudinal groove
thumb rest and catch
overall length 135 mm
straight, cross serrated jaws
curved, cross serrated jaws

G-17610
G-17615



ARRUGA
NEEDLE HOLDER
 14 mm heavy jaws
 thumb rest and catch
 overall length 140 mm
 straight, smooth jaws
 curved, smooth jaws

G-17640
 G-17645

G-17660 *HEIDELBERG MODEL*
NEEDLE HOLDER
 straight
 14 mm heavy jaws
 with ratchet lock
 overall length 125 mm



G-17650 *BOYNTON*
NEEDLE HOLDER
 straight
 with ratchet lock
 overall length 120 mm



G-17665 *ROHRSCHEIDER*
NEEDLE HOLDER
 straight
 13 mm heavy jaws
 1.3 x 1.5 mm tips
 overall length 100 mm



G-19270 *BARRAQUER*
RAZOR BLADE HOLDER
flat handle
overall length 110 mm



G-19320 *CASTROVIEJO*
RAZOR BLADE HOLDER
6 mm round handle
overall length 120 mm



G-19330 *TROUTMAN*
RAZOR BLADE HOLDER
straight
6 mm round handle
overall length 95 mm



G-19300 **RAZOR BLADES**
10 pcs. per box



G-18150 **HOLDER FOR ABSORBENT STRIPS**
serrated
to secure sterile absorbent swabs (G-18127)
overall length 175 mm

