

SCLERAL DEPRESSORS LOCALIZERS





G-32152 *LAQUA*
SCLERAL DEPRESSOR
with pocket clip
3 x 1.8 mm teardrop-shaped tip
overall length 128 mm



G-32150 *LAQUA*
SCLERAL DEPRESSOR
double-ended with pocket clip
3 x 1.8 mm teardrop-shaped tip
2 mm diameter scleral marker
overall length 134 mm



G-32153 *LUCKE*
SCLERAL DEPRESSOR
with pocket clip
3.8 x 2.0 mm teardrop-shaped tip
overall length 128 mm



G-32027 *KREISSIG*
SCLERAL DEPRESSOR
double-ended with pocket clip
conical tip with 2.4 mm diameter
and 6.0 mm bar
overall length 128 mm



G-15710 *MEYER-SCHWICKERATH*
LOCALIZER AND SCLERAL DEPRESSOR
 for fixation and rotation




G-32715 *SCHOCKET*
SCLERAL DEPRESSOR
 double-ended with pocket clip
 3 x 2 mm teardrop tip and
 6.0 mm bar
 overall length 132 mm



G-32714
SCLERAL DEPRESSOR
 double-ended
 2 mm and 4 mm diameter
 diamond dusted tips
 overall length 120 mm

SCLERAL DEPRESSORS · LOCALIZERS

13



G-37445 *HATTENBACH*
TRIANGLE SCLERAL DEPRESSOR
 for gentle buckling and firm grip without the risk of rupture
 overall length 133 mm