

PROBES DILATORS



G-150301 *WILDER*
LACRIMAL PROBE
delicate, long taper
0.45 mm tip diameter
overall length 104 mm



G-150302 *WILDER*
LACRIMAL PROBE
medium, medium taper
0.50 mm tip diameter
overall length 104 mm



G-150303 *WILDER*
LACRIMAL PROBE
strong, short taper
0.55 mm tip diameter
overall length 104 mm



G-150401 *WILDER*
LACRIMAL PROBE
slim
delicate, long taper
0.30 mm tip diameter
overall length 95 mm



G-150402 *WILDER*
LACRIMAL PROBE
slim
medium, medium taper
0.35 mm tip diameter
overall length 95 mm



G-150403 *WILDER*
LACRIMAL PROBE
slim
strong, short taper
0.40 mm tip diameter
overall length 95 mm



G-15020 *GIESSEN MODEL*
LACRIMAL PROBE
11 mm tapered shaft, angled
0.4 mm tip width
overall length 77 mm



KELLNAR
RETROGRADE PROBE
 for canalicular repair
 0.5 mm suture hole
 overall length 95 mm
 right
 left


G-15050
 G-15055



G-15095
PIGTAIL PROBE
 double-ended
 for canalicular repair
 with 0.4 mm suture holes
 overall length 140 mm




G-15096
PIGTAIL PROBE
 double-ended
 for canalicular repair
 with 1.0 mm suture hooks
 overall length 140 mm



PROBE FOR STROMAL CHANNEL
 for corneal implants in patients with
 keratoconus
 overall length 115 mm
 bent to the left
 bent to the right

G-37602
 G-37604





BOWMAN
LACRIMAL PROBE
cylindrical, silver
overall length 130 mm

- G-150001** size 0000/000
- G-150002** size 00/0
- G-150003** size 1/2
- G-150004** size 3/4
- G-150005** size 5/6



BOWMAN
LACRIMAL PROBE
cylindrical, stainless steel
overall length 130 mm

- G-150101** size 0000/000
- G-150102** size 00/0
- G-150103** size 1/2
- G-150104** size 3/4
- G-150105** size 5/6



BOWMAN
LACRIMAL PROBE
bulbous tips, silver
overall length 130 mm

- G-150051** size 0000/000
- G-150052** size 00/0
- G-150053** size 1/2
- G-150054** size 3/4
- G-150055** size 5/6



BOWMAN
LACRIMAL PROBE
bulbous tips, stainless steel
overall length 130 mm


- G-150151** size 0000/000
- G-150152** size 00/0
- G-150153** size 1/2
- G-150154** size 3/4
- G-150155** size 5/6

JUENEMANN

LACRIMAL PROBE

with 72 mm probe and 18 mm grip plate to insert suture material provided with cleaning adaptor G-31022

- G-33193 0.6 mm diameter
- G-33194 0.8 mm diameter
- G-33195 1.0 mm diameter



G-15315

KRUMEICH

DALK PROBE

for deep lamellar keratoplasty
12 mm tapered probe
overall length 73 mm



MACKENSEN

TRABECULAR PROBE

7 mm curved probe
overall length 92 mm

- G-15105 right
- G-15110 left



HARMS

TRABECULAR PROBE

dual 9.5 mm curved probes
3 mm spread
overall length 46 mm

- G-15125 right
- G-15130 left



14

THOMAS NEUHANN

U-SHAPED PROBE

with holding plate
dual 9 mm curved probes
3 mm spread
overall length 12 mm

- G-15136 0.30 mm probe diameter
- G-15138 0.25 mm probe diameter



G-15137

THOMAS NEUHANN

HOLDER

with lock
for U-shaped probes G-15136 and G-15138
overall length 92 mm

