



CANNULAS I/A INSTRUMENTS





LACRIMAL.....	189
ANTERIOR CHAMBER MAINTAINERS.....	191
HYDRODISSECTION.....	192
IRRIGATION VECTIS.....	193
VISCOELASTIC.....	194
CAPSULE POLISHERS.....	195
LASIK.....	196
BIMANUAL IRRIGATION/ASPIRATION.....	198
COAXIAL IRRIGATION/ASPIRATION.....	203
BACKFLUSH.....	206
INFUSION.....	209


LACRIMAL CANNULA
conical, straight
tapers from 19 to 25 gauge
25 mm shaft
silver
G-33140
G-33135 stainless steel



LACRIMAL CANNULA
conical, strongly curved
tapers from 19 to 25 gauge
25 mm shaft
silver
G-15140
G-15141 stainless steel



LACRIMAL CANNULA
conical, angled at 10 mm
tapers from 19 to 25 gauge
25 mm shaft
silver
G-15150
G-15151 stainless steel



LACRIMAL CANNULA
conical, angled at 18 mm
tapers from 19 to 25 gauge
25 mm shaft
silver
G-15160
G-15161 stainless steel


G-15190 *ANEL*
LACRIMAL CANNULA
straight
15 mm shaft
23 gauge / 0,64 mm

G-15195 *ANEL*
LACRIMAL CANNULA
strongly curved
16 mm shaft
23 gauge / 0,64 mm


LACRIMAL CANNULA
 cylindrical, straight
G-33160 20 gauge / 0.9 x 35 mm
G-33170 no. 12, 0.72 x 32 mm
G-33175 no. 14, 0.67 x 30 mm
G-33180 no. 16, 0.60 x 22 mm
G-33185 no. 18, 0.45 x 22 mm
G-33190 no. 20, 0.42 x 22 mm



G-15170 *BONN MODEL*
LACRIMAL CANNULA
 with probe tip
 1.15 mm diameter olive tip




BANGERTER
LACRIMAL PROBE CANNULA
 straight, 0.6 mm side port
G-15271 23 gauge / 0.64 x 60 mm
G-15272 21 gauge / 0.80 x 40 mm
G-15273 23 gauge / 0.64 x 40 mm




G-33197 *JUENEMANN*
LACRIMAL PROBE CANNULA
 to insert suture material
 21 gauge / 0.8 x 35 mm




G-33198 *JUENEMANN*
LACRIMAL PROBE CANNULA
 with grip plate
 65 mm probe
 to insert suture material
 21 gauge / 0.8 x 65 mm




G-15375 **RETROBULBAR INJECTION CANNULA**
G-32605 24 gauge / 0.6 x 60 mm
G-32600 24 gauge / 0.6 x 35 mm
G-15376 25 gauge / 0.5 x 35 mm
G-15377 25 gauge / 0.5 x 50 mm
G-15377 25 gauge / 0.5 x 40 mm




G-15200 *CRANZ*
EYE DOUCHE
 irrigation through five front ports
 silver



G-33416 *BOPP*
ANTERIOR CHAMBER INFUSION CANNULA
 anterior chamber maintainer with
 threaded tip for better fixation
 and silicone tube
 20 gauge / 0.9 mm




G-33422 *BONN MODEL*
ANTERIOR CHAMBER INFUSION CANNULA
 self-retaining textured tip
 dual 0.35 mm side ports
 silicone tube 15 cm
 23 gauge / 0.6 mm



G-16105 **IRRIGATION CYSTOTOME CANNULA**
 26 gauge / 0.45 mm




G-15220 **ENZYME CANNULA**
 24 gauge / 0.6 mm



RETROBULBAR · EYE DOUCHE · A.C. MAINTAINERS · CYSTOTOME · ENZYME


G-15275 CYCLODIALYSIS SPATULA CANNULA

12 mm long spatula
anterior 0.2 mm port
26 gauge / 0.45 mm



G-31355 BUSSE CYCLODIALYSIS SPATULA CANNULA



12 mm long spatula
dual 0.2 mm port
26 gauge / 0.45 mm



SAUTTER HYDRODISSECTION CANNULA

front opening
flattened tip
vertical, 27 gauge / 0.4 mm
vertical, 26 gauge / 0.45 mm
horizontal, 27 gauge / 0.4 mm

G-15250
G-32166
G-32167



vertical
G-15250
G-32166

horizontal
G-32167



G-32155 TUEBINGEN MODEL HYDRODISSECTION CANNULA

30 gauge / 0.3 mm




G-15185 FRANKFURT MODEL HYDRODISSECTION CANNULA

27 gauge / 0.4 mm



G-31895 HYDRODISSECTION CANNULA

vertically flattened tip
front opening
24 gauge / 0.6 x 20 mm



SPIRIG HYDRODISSECTION CANNULA


for manual lens fragmentation
right
left

G-15252
G-15253




G-16115 IRRIGATION VECTIS CANNULA

4 front ports
5 mm wide lens loop
curved




G-32759 IRRIGATION VECTIS CANNULA

small model
4 front ports
4 mm wide lens loop
curved




G-31510 KLÖCKNER EXPRESSOR CANNULA

8 mm wide lens loop
4 front ports
lightly curved




G-31807 HEALON - VISCOELASTIC CANNULA

flattened tip, angled at 6 mm
30 mm shaft
21 gauge / 0.8 mm




G-22619 PH. JACOBI TRABECULAR ASPIRATION CANNULA

angled at 13 mm
24 gauge / 0,6 mm




G-15188 G-15189 BURK IRRIGATION TRABECULOTOMY CANNULA

8 mm curved irrigation probe
30 gauge / 0,3 mm
angled right
angled left



G-33199 NAGATA SUB-TENON INJECTION CANNULA

for atraumatic penetration and
injection of anesthetic agent
24 gauge / 0.6 mm




IRRIGATION VECTIS · HEALON · TRABECULOTOMY · SUB-TENON


G-15180 **IRRIGATION CANNULA**
 bulbous tip, 0.75 mm diameter tip
 extra delicate
 26 gauge / 0.45 mm




G-15230 *BONN MODEL*
IRRIGATION CANNULA
 flattened tip, front opening
 delicate, 0.5 x 1.5 mm
G-15240 medium, 0.6 x 1.9 mm




G-16134 *BINKHORST*
IRRIGATION CANNULA
 U-shaped
 22 gauge / 0.7 mm
G-16135 straight
G-16136 angled right
 angled left




G-32760 *MUELLER*
IRRIGATION CANNULA
 J-shaped with probe tip
 26 gauge / 0.45 x 19 mm
G-32765 angled left
 angled right



G-16138 *RENTSCH*
IRRIGATION CANNULA
 J-shaped
 angled at 10 mm
 26 gauge / 0.45 x 24 mm




G-16106 *SIMCOE*
IRRIGATION / ASPIRATION CANNULA
 dual shaft
 0.3 mm diameter aspiration port
 front opening irrigation port
 23 gauge / 0.64 mm
 with silicone tube and Luer-Lock hub




G-16120 **IRRIGATION / ASPIRATION CANNULA
COAXIAL**

straight
aspiration port 0.4 mm
infusion tube 18 gauge / 1.2 mm




G-15225 **ASPIRATION CANNULA**
lightly curved, front opening
G-15226 standard, 17 gauge / 1.5 mm




G-16125 **KRATZ**
CAPSULE POLISHING CANNULA
angled at 2 mm
half-textured tip
23 gauge / 0.64 mm




G-16133 **TEPING**
ASPIRATION POLISHING CANNULA
0.3 mm posterior port
half-textured tip, curved
23 gauge / 0.64 mm




G-16130 **CAPSULE POLISHING CANNULA**
with half-textured ball tip
18 gauge / 1.2 mm



G-16132 **PAPST**
CAPSULE POLISHING CANNULA
0.75 mm hemisphere textured tip
25 gauge / 0.5 mm




G-32529 **KYU**
CAPSULE POLISHING CANNULA
2 mm disc-shaped tip
textured
G-32527 22 gauge / 0.7 mm
23 gauge / 0.64 mm




G-32528

GREITE
CAPSULE POLISHING CANNULA
 curette tip
 textured
 22 gauge / 0.7 mm



G-32774

BRAUWEILER
ASPIRATION CANNULA
 angled at 13 mm
 textured polisher tip
 0.3 mm beveled port
 21 gauge / 0.8 mm




G-32769

H.-R. KOCH
ASPIRATION CANNULA
 angled at 14 mm
 0.3 mm anterior port
 half-textured, spatulated tip
 21 gauge / 0.8 mm




G-32776

BRAUWEILER
IRRIGATION CANNULA
 angled at 12 mm
 front opening
 21 gauge / 0.8 mm




G-32772

BRAUWEILER
IRRIGATION CANNULA
 angled at 13 mm
 textured 0.8 mm iris retractor and
 capsule polisher
 21 gauge / 0.8 mm




G-33909

W. BAUER
LASIK CANNULA
 double-armed
 6 ports
 for irrigation of the flap and stromal bed
 25 gauge / 0.5 mm



G-33910

ZIRM
LASIK DISSECTOR CANNULA
flattened
front opening
25 gauge / 0.5 mm




G-33911

ZIRM
LASIK IRRIGATION CANNULA
flattened and textured tip
with 5 irrigation ports
20 gauge / 0.9 mm



G-15310

KRUMEICH
CORNEAL IRRIGATION CANNULA
for big bubble technique during DALK
semi-sharp dissector tip
0.2 mm posterior port
27 gauge / 0.4 mm



KRUMEICH
CORNEAL IRRIGATION CANNULA
for big bubble technique during DALK
angled at 12 mm
23 gauge / 0.6 mm
G-15312 0.4 mm anterior port
G-15314 0.4 mm posterior port



G-22550

IRRIGATION CANNULA HANDPIECE
overall length 115 mm



G-32015

IRRIGATION CANNULA HANDPIECE
male / female connectors
overall length 112 mm



G-32771 *BRAUWEILER*
IRRIGATION CANNULA HANDPIECE
 male / female connectors
 78 mm



G-32773 *BRAUWEILER*
ASPIRATION CANNULA HANDPIECE
 male / male connectors
 80 mm




G-32778 *KEERL*
IRRIGATION CANNULA HANDPIECE
 with silicone tube



G-32010 **IRRIGATION CANNULA HANDPIECE**
 with silicone tube



SCHMACK
 21 gauge / 0.8 mm




G-22632 **IRRIGATION HANDPIECE**
 0.5 mm port

G-22633 **ASPIRATION HANDPIECE**
 0.3 mm port




D.D. KOCH - T. KOHNEN
ASPIRATION HANDPIECE
 0.3 mm aspiration port
 textured tip



G-22100 22 gauge / 0.7 mm
G-22102 23 gauge / 0.6 mm


G-22101

D.D. KOCH - T. KOHNEN
IRRIGATION HANDPIECE
front opening
22 gauge / 0.7 mm




G-22105
G-22106

ASPIRATION HANDPIECE
0.25 mm aspiration port
conical tip
oval tube 0.95 x 0.55 mm
smooth
textured



G-22107

IRRIGATION HANDPIECE
0.5 mm irrigation port
conical tip
oval tube 0.95 x 0.55 mm




G-22637

SCHMACK
ASPIRATION HANDPIECE
4 x 0.15 mm posterior ports
textured tip
23 gauge / 0.6 mm




G-32011 ASPIRATION HANDPIECE

0.3 mm aspiration port
textured, conical tip
21 gauge / 0.8 mm




G-32012 IRRIGATION HANDPIECE

dual 0.5 mm irrigation ports
textured, conical tip
21 gauge / 0.8 mm



G-32018 BRAUWEILER ASPIRATION HANDPIECE

angled at 13 mm
0.3 mm beveled port
textured tip
21 gauge / 0.8 mm



G-32019 H.-R. KOCH ASPIRATION HANDPIECE

angled at 14 mm
0.3 mm anterior port
half-textured, spatulated tip
21 gauge / 0.8 mm



G-32016

BRAUWEILER
IRRIGATION HANDPIECE
angled at 12 mm
front opening
21 gauge / 0.8 mm



G-32017

BRAUWEILER
IRRIGATION HANDPIECE
angled at 13 mm
textured 0.8 mm iris retractor and
capsule polisher
21 gauge / 0.8 mm



G-32026

R. KOCH
ASPIRATION HANDPIECE
0.3 mm aspiration port
textured tip
19 gauge / 1.1 mm



G-32028

R. KOCH
IRRIGATION HANDPIECE
angled
dual 0.75 mm side ports
19 gauge / 1.1 mm




G-32014

R. KOCH

IRRIGATION CHOPPER

angled
dual 0.75 mm side ports
1 mm angled chopper tip
19 gauge / 1.1 mm



G-32059

DICK

IRRIGATION CHOPPER

angled
dual 0.75 mm side ports
1 mm curved chopper tip
19 gauge / 1.1 mm




G-38601

ALTHAUS / CARTSBURG

IRRIGATION DESCMET HOOK

hook diameter: 0.25 mm
hook length: 0.6 mm
for Descemetorhexis
rhombic 0.9 x 0.6 mm irrigation tube
overall length 98 mm




G-38602

THOMAS NEUHANN

IRRIGATION DESCMET SCRAPER

scraper width: 0.8 mm
scraper length: 0.6 mm
rhombic 0.9 x 0.6 mm irrigation tube
for peripheral scraping of the recipient bed
overall length 94 mm




G-38631

SZURMAN

DMEK IRRIGATION HANDPIECE

curved, closed tip
side port 0.3 mm
tube 23 gauge / 0.6 mm




G-38632

SZURMAN



DMEK IRRIGATION HANDPIECE

curved, front opening
tube 23 gauge / 0.6 mm



PEP 2.2 MM I/A HANDPIECE

Pure Efficiency Phaco
for 2.2 mm incisions
0.3 mm aspiration port
coaxial I/A straight tip
coaxial I/A 30° angled tip



G-24318


G-24418

with G-24311 PEP infusion sleeve

G-24311

PEP 2.2 MM INFUSION SLEEVE

Pure Efficiency Phaco
for 2.2 mm incisions
purple translucent silicone




COAXIAL IRRIGATION/ASPIRATION · DMEK IRRIGATION



PEP 2.4 MM I/A HANDPIECE

Pure Efficiency Phaco
for 2.4 mm incisions
0.3 mm aspiration port
coaxial I/A straight tip
coaxial I/A 30° angled tip



G-24332


G-24432

with G-24314 PEP infusion sleeve

G-24314

PEP 2.4 MM INFUSION SLEEVE

Pure Efficiency Phaco
for 2.4 mm incisions
green translucent silicone




COAXIAL IRRIGATION/ASPIRATION



G-24330

I/A HANDPIECE

coaxial straight tip
0.3 mm aspiration port
with G-24100 infusion sleeve




G-24430

G-24435

G-24440


I/A HANDPIECE

coaxial 30° angled tip
0.3 mm aspiration port
0.35 mm aspiration port
0.4 mm aspiration port





with G-24100 infusion sleeve


G-24100 **INFUSION SLEEVE**
blue silicone
for I/A handpieces




PAPST
I/A HANDPIECE
coaxial
infusion tube 16 gauge / 1.6 mm
aspiration port 0.3 mm, straight
G-22525
G-22527 aspiration port 0.4 mm, straight
G-22540 aspiration port 0.3 mm, angled 30°
G-22542 aspiration port 0.4 mm, angled 30°
G-22545 aspiration port 0.5 mm, angled 30°



G-22615 *PAPST*
I/A HANDPIECE
coaxial 30° angled tip
aspiration port 0.3 mm
infusion tube 18 gauge / 1.3 mm



G-22556 *HUNOLD*
I/A HANDPIECE
coaxial, curved
aspiration port 0.3 mm
infusion tube 18 gauge / 1.3 mm



COAXIAL IRRIGATION/ASPIRATION

G-33045 **BACKFLUSH HANDPIECE**
with Luer cone
overall length 110 mm



G-31805 **BACKFLUSH HANDPIECE**
with silicone chamber and Luer cone
overall length 110 mm



G-37002 **BACKFLUSH HANDPIECE**
with silicone chamber and
Luer-Lock connector
overall length 115 mm



G-37000 **BACKFLUSH HANDPIECE**
with silicone chamber and chuck
overall length 110 mm



for backflush tips G-37010 up to G-37080

SILICONE CHAMBER
for backflush handpieces
G-31805, G-37000, G-37002



G-31806 reusable
G-34806 single use, 10 pcs. per box, sterile

RETINAL · BACKFLUSH HANDPIECES

15

G-37005 **BACKFLUSH HANDPIECE**
with large silicone chamber
rotatable chuck



for backflush cannulas G-37010 up to G-37080



G-37006 **SILICONE CHAMBER**
large, for backflush handpiece G-37005





G-33050 **BACKFLUSH CANNULA**
35 mm shaft
20 gauge / 0.9 mm

G-33054
23 gauge / 0.6 mm

G-33053
25 gauge / 0.5 mm



G-33070 **BACKFLUSH CANNULA**
35 mm shaft
20 gauge / 0.9 mm shaft
25 gauge / 0.5 mm tip
3.5 mm tapered tip



G-33071 **BACKFLUSH CANNULA**
35 mm shaft
0.5 mm side port
20 gauge / 0.9 mm






G-33072 *LAQUA*
BACKFLUSH CANNULA
35 mm shaft
inside diameter: 21 gauge / 0.8 mm
outside diameter: 19 gauge / 1.1 mm

G-33073 *KAMPIK*
BACKFLUSH CANNULA
35 mm shaft
inside diameter: 22 gauge / 0.7 mm
outside diameter: 20 gauge / 0.9 mm


















G-33075 **BACKFLUSH CANNULA**
35 mm shaft
0.8 mm silicone brush
20 gauge / 0.9 mm





G-33077 **BACKFLUSH CANNULA**
20 gauge / 0.9 mm shaft
24 gauge / 0.55 mm tip
7 mm tapered tip, angled at 3.5 mm





<p>G-37010 G-37011 G-37012</p>	<p>BACKFLUSH TIP 35 mm shaft 20 gauge / 0.9 mm 23 gauge / 0.6 mm 25 gauge / 0.5 mm</p>		
<hr/>			
<p>G-37015</p>	<p>BACKFLUSH TIP 35 mm shaft angled at 3 mm 20 gauge / 0.9 mm</p>		
<hr/>			
<p>G-37020</p>	<p>BACKFLUSH TIP 35 mm shaft 0.8 mm silicone brush 20 gauge / 0.9 mm</p>		
<hr/>			
<p>G-37025</p>	<p>BACKFLUSH TIP 35 mm shaft angled at 6 mm 0.8 mm silicone brush 20 gauge / 0.9 mm</p>		
<hr/>			
<p>G-37035</p>	<p>BACKFLUSH TIP 35 mm shaft 0.8 mm silicone tip 20 gauge / 0.9 mm</p>		
<hr/>			
<p>G-37046</p>	<p><i>KLEIN</i> BACKFLUSH TIP 20 gauge / 0.9 mm shaft 26 gauge / 0.45 mm tip 2 mm tapered tip overall length 35 mm</p>		
<hr/>			
<p>G-37065</p>	<p>BACKFLUSH TIP 35 mm shaft 20 gauge / 0.9 mm shaft 30 gauge / 0.3 mm tip tapered tip, angled at 5 mm oval aspiration port</p>		


G-37080 **BACKFLUSH TIP**
35 mm shaft
20 gauge / 0.9 mm shaft
23 gauge / 0.6 mm tip
2 mm tapered tip
overall length 35 mm




INFUSION CANNULA
silicone tube 40 cm
Luer-Lock metal adapter
20 gauge / 0.9 mm
G-33420 beveled, 4.0 mm (single-use G-33463)
G-33425 beveled, 5.0 mm
G-33430 beveled, 6.0 mm




INFUSION CANNULA
silicone tube 40 cm
Luer-Lock plastic adapter
20 gauge / 0.9 mm
G-33450 beveled, 4.0 mm
G-33451 beveled, 5.0 mm
G-33452 beveled, 6.0 mm




INJECTION CANNULA
for viscous fluids, without silicone tube
for butterfly-19 cannula
20 gauge / 0.9 mm
G-33484 beveled, 4.0 mm
G-33485 beveled, 5.0 mm
G-33486 beveled, 6.0 mm




INJECTION CANNULA
for viscous fluids, without silicone tube
for butterfly-21 cannula
20 gauge / 0.9 mm
G-33494 beveled, 4.0 mm
G-33495 beveled, 5.0 mm
G-33496 beveled, 6.0 mm



DUISBURG MODEL
INJECTION CANNULA
for viscous fluids
25 cm silicone tube with 2 metal sleeves
and Luer-Lock adapter
20 gauge / 0.9 mm
G-32180 beveled, 4.0 mm
G-32181 beveled, 5.0 mm
G-32182 beveled, 6.0 mm



HAMBURG MODEL
INJECTION CANNULA
 for viscous fluids
 25 cm silicone tube with metal sleeve
 and Luer-Lock adapter
 3 spare silicone tubes




- G-33470** 20 gauge / 0.9 mm, beveled, 4 mm
- G-33471** 20 gauge / 0.9 mm, beveled, 5 mm
- G-33472** 20 gauge / 0.9 mm, beveled, 6 mm
- G-33473** 23 gauge / 0.6 mm, beveled, 4 mm
- G-33474** 23 gauge / 0.6 mm, beveled, 6 mm

G-33475 **SPARE SILICONE TUBE**
 for injection cannulas
 25 cm
 3 pcs. per box




WEBER
ASPIRATION CANNULA
 for viscous fluids
 inside diameter: 22 gauge / 0.7 mm
 outside diameter: 20 gauge / 0.9 mm




- G-33051** 5 mm shaft
- G-33052** 10 mm shaft


G-37390 **ASPIRATION CANNULA**
 23 gauge / 0.6 mm
 10 mm shaft




G-33056 *ROIDER*
ASPIRATION CANNULA
 for viscous fluids
 0.7 mm side port
 19 gauge / 1.0 mm



G-33055 **PFCL CANNULA**
 with dual shaft
 24 gauge / 0.55 mm
 25 gauge / 0.5 mm, lateral tube



G-33057 *CHANG*
PFCL CANNULA
 for injection of heavy fluids
 dual bore, coaxial
 25 gauge / 0.5 mm tip
 20 gauge / 0.9 mm shaft



HEIDELBERG MODEL

CANNULA


for injection or aspiration of viscous fluids
and Densiron[®] 68, bevel 30°

G-32699

19 gauge / 1.1 mm

G-32698

18 gauge / 1.2 mm



STOPPER

for viscous fluid aspiration
with tube connection for single-use syringe

G-33065

10 ml syringe

G-33066

20 ml syringe



PRESSURE TUBE


for injection of viscous fluids
Luer-Lock female / male

G-32697

reusable

G-32696

single use, 10 pcs. per box, sterile



G-28752

SYRINGE HOLDER

for manual injection of
viscous fluids in glass syringes,
with clamp and retraction mechanism

SYRINGE CLIP
customized labeling
for 2 ml syringes
for 5 ml syringes

G-19960
G-19961



SILICONE BULB
with Luer Cone
20 ml
60 ml

G-32665
G-32670

