


MISCELLANEOUS




CAPSULE TENSION RING INJECTORS	241
PUPIL DILATOR INJECTOR	241
LACRIMAL / DCR	242
LASEK	243
DIAMOND BURR	244
VITRECTOMY LENSES	244
EYE SHIELD	245
EYE SPHERES	245
CAUTERY	246




G-32955 **INJECTOR**
for capsular tension ring
shaft diameter 1.2 mm
for 1.9 mm incisions
shaft length 30 mm
overall length 156 mm




G-32960 **INJECTOR**
for *MORCHER* capsular tension ring
shaft diameter 1.2 mm
for 1.9 mm incisions
shaft length 30 mm
overall length 156 mm



G-32965 **INJECTOR**
for *MORCHER* capsular ring
shaft diameter 1.6 mm
for 2.4 - 2.5 mm incisions
shaft length 30 mm
overall length 156 mm



G-32970 **INJECTOR**
for *MORCHER* pupil dilator
overall length 140 mm



G-31810 *FRIEDMANN*
RONGEUR
overall length 140 mm



G-32355 **NASAL SPECULUM**
overall length 150 mm



G-33125 *MEAD*
HAMMER
300 g, lead filled
overall length 190 mm



G-18630 *HEIDELBERG MODEL*
ELEVATORIUM
overall length 143 mm



G-18620 *HEIDELBERG MODEL*
RASPATORY
overall length 155 mm






G-18240 *WEST*
BONE CHISEL
for lacrimal sac
flat, 5.0 mm wide
overall length 160 mm




G-18245 *WEST*
BONE GOUGE
for lacrimal sac
concave, 5.0 mm wide
overall length 162 mm




ELSA
LASEK TREPHINE
with crosshairs
overall length 100 mm
G-33650 8.0 mm
G-33652 9.0 mm



LOHMANN
LASEK FUNNEL
with crosshairs
overall length 100 mm
G-33655 8.5 mm
G-33657 9.5 mm



G-32313 *WRIGHT*
FASCIA NEEDLE
overall length 115 mm



DIAMOND BURR

top, inside, outside coat

- G-18300 8.0 mm
- G-18301 9.0 mm
- G-18302 10.0 mm



- G-36000 *KILP*
LENS
60 D flat, field of vision 20°, silicone ring
central posterior pole and central vitreous body



- G-36010 *LANDERS*
LENS
93 D, biconcave with silicone ring
clear view of fundus in air-filled phakic eye



- G-36015 *MACHEMER*
LENS
with silicone ring to magnify structures
in vitreous cavity or on retinal surface



- G-36020 *PEYMAN*
LENS
wide field of vision 30°, with silicone ring
for central fundus and central vitreous body



- G-36025 *TOLENTINO*
LENS
prism 20° with silicone ring
posterior peripheral fundus and
posterior peripheral vitreous body



- G-36030 *TOLENTINO*
LENS
prism 30° with silicone ring
peripheral viewing beyond the equator



- G-36045 **SILICONE REPLACEMENT RING**
for *KILP* lens G-36000



G-36050 **SILICONE REPLACEMENT RING**
for lenses from G-36010 to G-36035



G-36065 **LENS HANDLE RING**
for lenses without silicone ring
to suture on to sclera
for G-36000



G-31105 **EYESHIELD**
plastic, transparent for use on both eyes
100 pcs. per bag




CIBIS
EYE SPHERE
G-20216 16 mm
G-20217 17 mm
G-20218 18 mm
G-20219 19 mm
G-20220 20 mm



HARMS
EYE SPHERE
G-20316 16 mm
G-20317 17 mm
G-20318 18 mm
G-20319 19 mm
G-20320 20 mm



WALSER
EYE SPHERE
G-20416 16 mm
G-20417 17 mm
G-20418 18 mm
G-20419 19 mm
G-20420 20 mm



17
EYE SHIELD · EYE SPHERES

G-19820 **BATTERY CAUTERY**
with spare tip
and batteries
overall length with
spare tip 155 mm



9/10

G-19825 **SPARE TIP**
for battery cauterization G-19820
overall length 50 mm



G-31900 *KUHN*
INTRAOCULAR MAGNET
30 mm shaft, 1.9 mm outside diameter
4 mm magnetic tip
overall length 90 mm

