

STERILIZATION CONTAINERS AND ADAPTERS



MICROSTOP®	249
STEEL CASES	252
TRAYS & RACKS FOR STEEL CASES	256
SILICONE STRIPS & MATS	257
SMALL CASES	258
STERISAFE® DURO	259
ADAPTERS	264

G-39745
**MICROSTOP®
 CONTAINER LOW**

for one tray
 outside dimensions: 310 x 189 x 70 mm
 inside dimensions: 283 x 177 x 50 mm

**G-39750**
**MICROSTOP®
 CONTAINER HIGH**

for two trays
 outside dimensions: 310 x 189 x 90 mm
 inside dimensions: 283 x 177 x 70 mm


**MICROSTOP®
 TRAY MINI**

29 mm, low
 size: 277 x 172 x 29 mm
 without lid
 G-39783
 with lid
 G-39786


**MICROSTOP®
 TRAY MINI**
G-39785
G-39784

38 mm, high
 size: 277 x 172 x 38 mm
 without lid
 with lid



MICROSTOP®

18

G-39787

**MICROSTOP®
TRAY LID MINI**

size: 277 x 172 mm
for G-39783 and G-39785



G-39782

G-39789

MICROSTOP®

RETAINER ELEMENT

for MICROSTOP® tray mini
nap Style, 120 mm
separable into 6 single elements
29 mm, low
38 mm, high



G-39788

MICROSTOP®

CLIP FOR RETAINER ELEMENTS

G-39782 and G-39789
10 pcs. per bag



G-39791

MICROSTOP®

ID LABEL

for tray, aluminium
customized labeling
silver-colored



G-39741

MICROSTOP®

LABEL

for lid, customized labeling



G-39746

MICROSTOP®

ID LABEL

for front, customized labeling




MICROSTOP®
LOGISTIC FRAME
 for lid to use with labels

G-39753 color: red
G-39754 color: blue
G-39755 color: green



G-39762 **MICROSTOP®**
SEAL
 with indicator
 1000 pcs. per box



G-39764 **MICROSTOP®**
TAG
 with indicator
 1000 pcs. per box



G-39680 **PROTOCOL CARDS**
 for ERGO-SAFE tray
 320 pcs. per box



G-19822

STERILIZATION CASE

- for steam sterilization
outside dimensions: 220 x 165 x 55 mm
- textile filter in lid and base
 - silicone seal
 - perforated insert without edge
 - perforated silicone tray G-32785
 - 2 labels with clamp



G-19823

STERILIZATION CASE

- for steam sterilization
outside dimensions: 220 x 165 x 75 mm
- textile filter in lid and base
 - silicone seal
 - perforated insert without edge
 - perforated silicone tray G-32785
 - 2 labels with clamp



G-28151

STERILIZATION CASE

- for steam sterilization
outside dimensions: 300 x 220 x 75 mm
- textile filter in lid and base
 - silicone seal
 - perforated insert without edge
 - 3 silicone burl strips
 - 2 silicone paddings in lid
 - 2 labels with clamp



G-28150

STERILIZATION CASE

- for steam sterilization
outside dimensions: 300 x 220 x 75 mm
- textile filter in lid and base
 - silicone seal in the lid
 - perforated insert with edge, subdivided, with perforated silicone plate and perforated silicone tray
 - silicone insert, small, with perforated silicone tray
 - 2 labels with clamp



G-19876**STERILIZATION CASE**

- for steam sterilization
 outside dimensions: 300 x 220 x 55 mm
- textile filter in lid and base
 - silicone seal
 - perforated insert without edge
 - 3 silicone burl strips
 - 2 silicone paddings in lid
 - 2 labels with clamp

**G-19856****STERILIZATION CASE**

- for steam sterilization
 outside dimensions: 300 x 220 x 55 mm
- paper filter with quick-clasp in lid and base
 - silicone seal in the lid
 - without insert and profile strips
 - 2 labels with clamp

**G-19843****STERILIZATION CASE**

- for steam sterilization
 outside dimensions: 300 x 220 x 75 mm
- textile filter in lid and base
 - silicone seal
 - perforated insert with edge
 - perforated silicone plate
 - perforated insert without edge
 - 3 silicone burl strips
 - 2 silicone paddings in lid
 - 2 labels with clamp



G-19842

STERILIZATION CASE

- for steam sterilization
outside dimensions: 300 x 220 x 75 mm
- textile filter in lid and base
 - silicone seal in the lid
 - perforated insert with edge and perforated silicone plate
 - perforated insert with edge perforated silicone plate and silicone instrument tray
 - 2 labels with clamp



G-19877

STERILIZATION CASE

- for steam sterilization
outside dimensions: 400 x 290 x 55 mm
- textile filter in lid and base
 - silicone seal
 - perforated insert without edge
 - 4 silicone burl strips
 - 3 silicone paddings in lid
 - 2 labels with clamp



G-19892

STERILIZATION CASE

- for steam sterilization
outside dimensions: 400 x 290 x 70 mm
- textile filter in lid and base
 - silicone seal in the lid
 - 2 perforated inserts without edge with 4 silicone nap strips each
 - 3 silicone profile strips at the base of the upper perforated insert
 - 3 silicone profile strips in the lid
 - 2 labels with clamp



G-19894

STERILIZATION CASE

- for steam sterilization
- outside dimensions: 400 x 290 x 70 mm
- textile filter in lid and base
- silicone seal in the lid
- perforated insert without edge with 4 silicone nap strips
- perforated insert with edge and silicone nap plate
- 3 silicone profile strips in the lid
- 2 labels with clamp



G-19890

STERILIZATION CASE

- for steam sterilization
- outside dimensions: 495 x 270 x 70 mm
- textile filter in lid and base
- silicone seal in the lid
- perforated insert with edge 4 silicone nap strips and lateral compartment (240 x 110 x 40 mm),
- 3 silicone profile strips in the lid
- 2 labels with clamp



PERMANENT TEXTILE FILTER

- for steam sterilization cases

G-31650

- for G-19822, G-19823

G-31651

- for G-19841, G-19842, G-19843, G-19876, G-28150, G-28151

G-31652


- 6 perforations, for G-19877, G-19892, G-19894

G-31656

- for G-19877

G-31661

- for G-19890



G-31649

INDICATOR

- for steam sterilization, slim
- 1000 pcs. per box



- G-32675** **SILICONE TRAY**
with bottom
outside dimensions: 160 x 85 x 20 mm




- G-32660** **STERILIZATION RACK**
for 12 cannulas, silicone
outside dimensions: 85 x 56 x 8 mm



- G-19870** **STERILIZATION RACK**
for 10 instruments with handle, silicone
outside dimensions: 160 x 85 x 20 mm



- G-32785** **SILICONE PERFORATED TRAY**
outside dimensions: 194 x 135 x 24 mm



G-39650**STERILIZATION CASE**


for steam sterilization, plastic
 for 6 vitreoretinal instruments
 outside dimensions: 210 x 130 x 26 mm

**G-32676****SILICONE NAP STRIP BLUE**

340 x 30 mm

**G-32677****SILICONE STRIP BLUE**

350 mm long

**SILICONE BURL MAT**

perforated

G-19866


240 x 230 mm

G-19865


300 x 150 mm

G-19867

470 x 245 mm

**G-19899****SILICONE PLATE**

perforated, 270 x 185 mm
 for hot air sterilization cases
 G-19841, G-19842, G-19843,
 G-19846, G-19847



G-31657 **SAFETY SEAL**
for steam sterilization cases
100 pcs. per box



G-31659 **SEAL**
for hot-air sterilization cases
100 pcs. per box



G-32695 **SILICONE CLOSURE**
17 cm long



G-33080 **STERILIZATION CASE**
for steam sterilization, plastic, for 5 instruments with handle
outside dimensions: 160 x 85 x 20 mm



G-33090 **STERILIZATION CASE**
for one instrument or diamond knife, plastic
outside dimensions: 160 x 60 x 20 mm



G-33100 **STERILIZATION CASE**
for one vitreoretinal instrument, plastic
outside dimensions: 195 x 60 x 25 mm



G-36090

STERILIZATION CASE

plastic, for 6 vitreoretinal lenses



G-39527

STERISAFE® DURO A4 PLUS

container for steam, gas (EO, FO) and plasma sterilization, without trays and strips
outside dimensions: 290 x 190 x 70 mm
inside dimensions: 265 x 165 x 60 mm



G-39526

STERISAFE® DURO A4 PLUS COMPLETE SET

container for steam, gas (EO, FO) and plasma sterilization, including:
1 x Sterisafe® wire basket A4-30
5 x Sterisafe® silicone retainer strips 230 mm
1 x Sterisafe® metal tray A4-30
outside dimensions: 290 x 190 x 70 mm
inside dimensions: 265 x 165 x 60 mm

G-39525

STERISAFE® DURO A4

container for steam, gas (EO, FO) and plasma sterilization
outside dimensions: 290 x 190 x 60 mm
inside dimensions: 265 x 165 x 49 mm



G-39528

STERISAFE® DURO A6

container for steam, gas (EO, FO) and plasma sterilization
outside dimensions: 190 x 110 x 60 mm
inside dimensions: 165 x 92 x 49 mm



STERISAFE® DURO

18

G-39543 STERISAFE® DURO WIRE BASKET A4-30

outside dimensions: 260 x 163 x 30 mm

inside dimensions: 249 x 157 x 23 mm



G-39541 STERISAFE® DURO METAL TRAY A4-15

outside dimensions: 260 x 163 x 15 mm

inside dimensions: 249 x 157 x 11 mm



G-39542 STERISAFE® DURO METAL TRAY A4-30

outside dimensions: 260 x 163 x 30 mm

inside dimensions: 249 x 157 x 26 mm



G-39549 STERISAFE® DURO WIRE BASKET A6-30

outside dimensions: 162 x 90 x 30 mm

inside dimensions: 156 x 80 x 23 mm



G-39548 STERISAFE® DURO METAL TRAY A6-15

outside dimensions: 162 x 90 x 15 mm

inside dimensions: 156 x 80 x 11 mm



G-39545 STERISAFE® SILICONE RETAINER STRIP 75
for Sterisafe® DURO A4 and A6
length: 75 mm



G-39546 STERISAFE® SILICONE RETAINER STRIP 145
for Sterisafe® DURO A4 and A6
length: 145 mm



G-39547 STERISAFE® SILICONE RETAINER STRIP 230
for Sterisafe® DURO A4 and A6
length: 230 mm



G-39501 STERISAFE® FILTER
A4 10 pcs. per box
G-39511 A6 10 pcs. per box



G-39530 STERISAFE® FILTER TAG EXPIRATION DATE
10 pcs. per box




G-39535 SEAL / LABEL
seal with color indicator and blank label
2000 pcs. per box

Steri.- Datum/Date	
Name	
	Dampf (braun)  Steam (brown)


G-39540 **SECURITY SEAL**
100 pcs. per box



G-39502 **STERISAFE® SILICONE BURL MAT** A4
G-39512 A6



G-39503 **STERISAFE® SILICONE GRID INSERT** A4
G-39513 A6




G-39504 **STERISAFE® METAL GRID INSERT** A4
G-39514 A6



G-39505

STERISAFE® METAL GRID INSERT

with 2 silicone burl strips




G-31670

MELAG METAL GRID INSERT

with 2 silicone nap strips

outside dimensions: 370 x 147 x 50 mm

for MELAG sterilizing case



G-31028 **CLEANING ADAPTER**
for ultrasonic handpieces (irrigation connector)
Luer female / male



G-31029 **CLEANING ADAPTER**
for ultrasonic handpieces (aspiration connector)
Luer female / female




G-31021 **CLEANING ADAPTER**
for ultrasonic tips



G-31013 **CLEANING ADAPTER**
for *HATTENBACH* 25 gauge hybrid instruments
(G-32944, G-32947 and G-32948)

G-22636 **ADAPTER**
for *PAPST* I/A handpiece
light model



G-31025 **ASPIRATION ADAPTER**
metal
for cleaning I/A handpieces
Luer-Lock female / female



G-31026 **ASPIRATION ADAPTER**
plastic
for cleaning I/A handpieces
Luer-Lock female / female



G-31016 **CLEANING ADAPTER**
 for diamond knives, small



G-31017 **CLEANING ADAPTER**
 for diamond knives, large



G-31022 **CLEANING ADAPTER**
 for G-37010 - G-37080



G-31031 **CLEANING ADAPTER**
 for aspiration handpiece G-37100 and G-37105
 Luer-Lock female / silicone tube
 2.1 x 1.0 mm



G-31018 **CLEANING ADAPTER**
 for PFCL cannula G-33055
 Luer-Lock female / silicone tube



G-31024 **CLEANING ADAPTER**
 for eye specula with aspiration



G-31023 **CLEANING ADAPTER**
 for eye specula with aspiration
 G-33141, G-33142 and G-33143



G-31027 **CLEANING ADAPTER**
 for *JUENEMANN* probes





G-31034 **CLEANING ADAPTER**
for foreign body instruments G-14530, G-14540, G-14550
Luer-Lock female



G-31014 **CLEANING ADAPTER**
for cannulas G-33057 and G-33059
Luer-Lock female



G-36604 **CLEANING ADAPTER**
for rotatable and interchangeable tips

PLEASE READ YOUR INSTRUMENT'S
CLEANING INSTRUCTIONS OR ASK
YOUR LOCAL SALES REPRESENTATIVE
FOR THE RIGHT CLEANING ADAPTER.