

FOREIGN BODY INSTRUMENTS





G-14550 **FOREIGN BODY INSTRUMENT**
 double-ended with sliding handle
 2.2 mm wide curved needle
 1 mm wide *PLANGE* spud
 overall length 95 mm



G-14540 **FOREIGN BODY INSTRUMENT**
 double ended with sliding handle
 2.2 mm wide curved needle
 1.4 mm wide gouge
 overall length 95 mm



G-14530 **FOREIGN BODY INSTRUMENT**
 double ended with sliding handle
 2.2 mm wide curved needle
 1.4 mm wide spud
 overall length 95 mm

G-14490 **FOREIGN BODY NEEDLE**
straight
4.5 mm long triangular tip
2.4 mm wide
overall length 120 mm



G-14500 **FOREIGN BODY NEEDLE**
curved
4.5 mm long triangular tip
2.4 mm wide
overall length 120 mm



G-14570 *PLANGE*
FOREIGN BODY SPUD
curved and rounded
1 mm wide spud
overall length 125 mm



G-14480 **FOREIGN BODY GOUGE**
straight
1.4 mm wide gouge
overall length 116 mm

