


GLAUCOMA PUNCHES RONGEURS



G-18210

CITELLI
BONE NIBBLING RONGEUR
3 mm x 4.5 mm bite
overall length 150 mm




G-18220

BEYER
BONE NIBBLING RONGEUR
2 mm x 3 mm bite
overall length 140 mm



G-18225

KERRISON
BONE NIBBLING RONGEUR
5 mm x 3.5 mm bite
overall length 160 mm



08

RONGEURS




G-18238

KELLY

DESCEMET'S MEMBRANE PUNCH

for Glaucoma
19 gauge / 1.0 mm
0.9 mm x 0.45 mm bite
overall length 122 mm



08

GLAUCOMA PUNCHES




G-18239

LUNTZ-DODICK

TRABECULECTOMY PUNCH

for Glaucoma
with 360° rotatable cutting head
1 mm x 0.5 mm bite
overall length 115 mm




G-18235

GREHN

MICRO TRABECULECTOMY PUNCH


for Glaucoma
20 gauge / 0.9 mm
0.7 mm x 0.3 mm bite
overall length 90 mm





G-18237

GASS
CORNEOSCLERAL PUNCH
for Glaucoma
with 360° rotatable cutting head
1,5 mm x 0,6 mm bite
overall length 85 mm




1.5 x 0.6 mm



G-18232

MOHR
GLAUCOMA PUNCH
1 mm circular bite
overall length 78 mm



1.0 mm

08

GLAUCOMA PUNCHES

