

## SCISSORS



CORNEAL .....	<b>81</b>
STITCH .....	<b>82</b>
CAPSULE .....	<b>84</b>
IRIS .....	<b>86</b>
IOL .....	<b>87</b>
TENOTOMY .....	<b>89</b>
STRABISMUS .....	<b>90</b>
MISCELLANEOUS · RING HANDLE .....	<b>91</b>
ENUCLEATION .....	<b>94</b>
VITREORETINAL .....	<b>95</b>



*TROUTMAN*

**CORNEAL SCISSORS**

7 mm blades  
9.5 mm from tip to pivot  
overall length 100 mm  
right, curved blades, blunt tips  
left, curved blades, blunt tips

**G-19710**  
**G-19715**



*CASTROVIEJO*

**CORNEAL SCISSORS**

lower blade 9.5 mm  
upper blade 9.0 mm  
overall length 105 mm  
right, curved blades, blunt tips  
left, curved blades, blunt tips

**G-19690**  
**G-19695**



*KATZIN*

**KERATOPLASTY SCISSORS**

6 mm blades  
8 mm from tip to pivot  
overall length 100 mm  
right, curved blades, blunt tips  
left, curved blades, blunt tips

**G-19720**  
**G-19725**



*SAUTTER*

**CORNEAL SCISSORS**

light curve, fine blunt tips  
10 mm blades  
14 mm from tip to pivot  
overall length 100 mm

**G-19718**



*CASTROVIEJO*

**CORNEAL SCISSORS**

strong curve, blunt tips  
12 mm blades  
15 mm from tip to pivot  
overall length 103 mm

**G-19780**



*CASTROVIEJO*

**CORNEAL SCISSORS**

angled to side, blunt tips  
12 mm blades  
overall length 110 mm

**G-19790**



**CORNEAL SCISSORS**

light curve  
8 mm blades  
overall length 105 mm  
sharp tips  
blunt tips

G-19600  
G-19605



**CORNEAL SCISSORS**

strong curve, blunt tips  
6 mm blades, overall length 105 mm  
11 mm blades, overall length 110 mm

G-19607  
G-19615



*MACKENSEN*  
**KERATOPLASTY SCISSORS**

strong curve, blunt tips  
overall length 85 mm

G-19727



**CORNEAL SECTION SCISSORS**

angled up, 8 mm blades  
blunt tips  
overall length 100 mm

G-19585



**MICRO STITCH SCISSORS**

straight, 8 mm blades  
sharp tips  
overall length 105 mm

G-19570



*FINE*  
**STITCH SCISSORS**

straight, 14 mm blades  
sharp tips  
overall length 110 mm

G-19580



**G-19590** *CASTROVIEJO*  
**ANTERIOR SYNECHIA SCISSORS**  
angled up  
9 mm vaulted blades, blunt tips  
overall length 105 mm



**G-19740** *CASTROVIEJO*  
**ANTERIOR SYNECHIA SCISSORS**  
angled up  
15 mm vaulted blades, blunt tips  
overall length 115 mm



**TRABECULUM SCISSORS**

sharp tips  
8.5 mm blades  
1.5 mm pivot width  
overall length 80 mm



**G-19745** straight  
**G-19747** curved

**CAPSULOTOMY SCISSORS**

angled up, sharp tips  
2 mm pivot width  
overall length 80 mm  
**G-19776** 9 mm blades  
**G-19777** 11 mm blades



**G-32100** **CAPSULOTOMY SCISSORS**  
angled, lightly curved tips  
13 mm blades  
2 mm pivot width  
overall length 80 mm



**G-32101** **CAPSULOTOMY SCISSORS**  
angled, strongly curved tips  
13 mm blades  
2 mm pivot width  
overall length 80 mm



**G-32105** **CAPSULOTOMY SCISSORS**  
lightly curved, sharp tips  
11 mm blades  
1.8 mm pivot width  
overall length 88 mm



**G-32023**      **CAPSULOTOMY SCISSORS**  
 curved, sharp tips  
 7 mm blades  
 1.8 mm pivot width  
 overall length 80 mm



**G-19620**      *TUEBINGEN MODEL*  
**CAPSULOTOMY SCISSORS**  
 straight, sharp tips  
 9 mm blades  
 1.8 mm pivot width  
 overall length 88 mm



**G-19630**      *TUEBINGEN MODEL*  
**CAPSULOTOMY SCISSORS**  
 curved, sharp tips  
 9 mm blades  
 1.8 mm pivot width  
 overall length 88 mm



**G-19760**      *VANNAS*  
**CAPSULOTOMY SCISSORS**  
 extra delicate, straight  
 6 mm blades, sharp tips  
 1.7 mm pivot width  
 overall length 80 mm



**G-19765**      *VANNAS*  
**CAPSULOTOMY SCISSORS**  
 extra delicate, curved  
 6 mm blades, sharp tips  
 1.7 mm pivot width  
 overall length 80 mm



**G-19770**      *VANNAS*  
**CAPSULOTOMY SCISSORS**  
 extra delicate, angled to side  
 7 mm blades, sharp tips  
 overall length 80 mm




09

CAPSULOTOMY · CAPSULE


**G-19775** *VANNAS*  
**CAPSULOTOMY SCISSORS**  
extra delicate, angled up  
7 mm blades, sharp tips  
overall length 80 mm



**G-32940** *H.-R. KOCH*  
**CAPSULE SCISSORS**  
horizontal, 45° angled  
2.5 mm blades, 360° rotatable head  
22 gauge / 0.7 mm straight shaft  
overall length 107 mm



**G-32942** *BRAUWEILER*  
**CAPSULE SCISSORS**  
vertical, 1.25 mm blades  
22 gauge / 0.7 mm curved shaft  
overall length 100 mm



**G-38246** *YAE*  
**CAPSULE SCISSORS**  
interchangeable, titanium, gold-colored  
3 mm straight micro blades  
23 gauge / 0.6 mm curved shaft  
for cutting synechia, iridotomies,  
and to rescue the rhexis during MICS  
overall length 42 mm  
for use with handle G-38240



**G-38240** *JAPAN MODEL*  
**HANDLE**  
for capsule scissors and forceps  
6 mm, titanium  
overall length 90 mm




09


CAPSULOTOMY · ROTATABLE · INTERCHANGEABLE




*OSHER*  
**CAPSULE SCISSORS**  
 360° rotatable and interchangeable tips  
 22 gauge / 0.7 mm  
 overall length 45 mm  
 to use with handpiece G-38215



*OSHER*  
**CAPSULE SCISSORS**  
 45° angled  
**G-38232**




*OSHER*  
**CAPSULE SCISSORS**  
 horizontal, 45° angled  
**G-38236**



*OSHER*  
**CAPSULE SCISSORS**  
 horizontal, 90° angled  
**G-38238**

*OSHER*  
**HANDPIECE FOR CAPSULE SCISSORS**  
 8 mm titanium handpiece, gold-colored  
**G-38215**




09

CAPSULE (INTERCHANGEABLE AND ROTATABLE) · IRIS

**G-32948**

*HATTENBACH*  
**HYBRID SCISSORS**  
 to cut sutures during iris and dislocated IOL fixation  
 1.25 mm blades  
 25 gauge / 0.5 mm angled shaft



**G-19540**  
**G-19545**  
**G-19550**

*NOYES*  
**IRIS SCISSORS**  
 straight  
 18 mm blades  
 overall length 120 mm  
 sharp tips  
 sharp-blunt  
 blunt tips



**G-19640**  
**G-19655**

*WECKER*  
**IRIS SCISSORS**  
 11 mm blades  
 overall length 110 mm  
 sharp tips  
 blunt tips





*WECKER*  
**IRIS SCISSORS**  
7 mm blades  
overall length 110 mm  
sharp tips  
blunt tips


**G-19660**  
**G-19675**




**G-19680** *BARRAQUER*  
**IRIS SCISSORS**  
7 mm blades  
blunt tips  
overall length 56 mm



**G-32945** *STEINMETZ*  
**IOL SCISSORS**  
6.5 mm upper blade  
serrated lower tip with groove  
overall length 92 mm  
For cutting acrylic and silicone foldable IOLs



**G-32946** *STEINMETZ*  
**IOL SCISSORS**  
for small incision technique  
6.5 mm upper blade  
serrated lower tip with groove  
overall length 92 mm  
For cutting acrylic and silicone foldable IOLs






*SNYDER-OSHER*

**IOL SCISSORS**

For cutting acrylic and silicone foldable IOLs  
 lower notched blade to create IOL traction  
 18 gauge / 1.2 mm shaft  
 overall length 150 mm  
 short model, 3.5 mm blades  
 long model, 4.5 mm blades



**G-32998**

**G-32996**



**G-32800**

*H.-R. KOCH*

**SILICONE LENS BI-SECTOR**

for dividing implanted silicone lenses  
 overall length without wire loop 130 mm  
 5 wire loops included (G-32811)

**WIRE LOOP**

for *H.-R. KOCH* silicone lens bi-sector  
 (G-32800)  
 overall length 65 mm (loop 12 mm)



**G-32810**

**G-32811**

1 piece


5 pieces



**COMBINED**

**NEEDLE HOLDER / SCISSORS**

curved 12 mm blades with  
 3 mm needle holder and tying platform  
 used to hold, place, and cut fine suture material  
 overall length 125 mm



**G-17440**

**G-17441**

without lock

with lock

09

IOL - SCISSORS WITH NEEDLE HOLDER

**G-19750** *WESTCOTT*  
**SCISSORS**  
blunt tips  
curved, standard blades  
overall length 120 mm



**G-19751** *WESTCOTT*  
**SCISSORS**  
pointed tips  
curved, standard blades  
overall length 120 mm



**G-32250** *WESTCOTT*  
**SCISSORS**  
blunt tips  
curved, delicate blades  
overall length 100 mm



*STEVENS*  
**TENOTOMY SCISSORS**  
short blades, blunt tips  
22 mm from tip to pivot  
overall length 105 mm  
**G-19470** straight  
**G-19475** curved



*STEVENS*  
**TENOTOMY SCISSORS**  
long blades, blunt tips  
28 mm from tip to pivot  
overall length 110 mm  
**G-19480** straight  
**G-19485** curved



G-32614

BONN MODEL

**TENOTOMY SCISSORS**

curved blades, blunt tips  
24 mm from tip to pivot  
overall length 90 mm



G-19430

G-19435

BONN MODEL

**STRABISMUS SCISSORS**

blunt tips  
24 mm from tip to pivot  
overall length 90 mm  
straight  
curved



G-19420

G-19425

**STRABISMUS SCISSORS**

long blade  
28 mm from tip to pivot  
overall length 105 mm  
straight  
curved



09

TENOTOMY · STRABISMUS · RING HANDLE

**STRABISMUS SCISSORS**

short blade  
28 mm from tip to pivot  
overall length 110 mm  
straight  
curved

**G-19460**  
**G-19465**



**EYE SCISSORS**

straight, sharp tips  
overall length:      tip to pivot:

**G-19350**  
**G-19355**  
**G-19360**  
**G-19365**

90 mm                  27 mm  
100 mm                29 mm  
110 mm                30 mm  
120 mm                36 mm



**EYE SCISSORS**

curved, sharp tips  
overall length:      tip to pivot:

**G-19370**  
**G-19375**  
**G-19380**  
**G-19385**

90 mm                  25 mm  
100 mm                28 mm  
110 mm                29 mm  
120 mm                34 mm



**EYE SCISSORS**

delicate, sharp tips  
31 mm from tip to pivot  
overall length 105 mm  
straight  
curved

G-19390  
G-19395



*BONN MODEL*  
**EYE SCISSORS**

extra delicate  
sharp tips  
19 mm from tip to pivot  
overall length 90 mm  
curved  
straight

G-19400  
G-19405



G-19410

*SPENCER*  
**STITCH SCISSORS**  
20 mm straight blades  
with hook to lift and  
secure sutures before cutting  
overall length 90 mm



09

UTILITY · STITCH · RING HANDLE

**G-19450** *GRAEFE*  
**EYE SCISSORS**  
angled up  
blunt tips  
overall length 90 mm



**G-19505** **EYE SCISSORS**  
angled to side, probe tip  
sharp-blunt  
overall length 115 mm



**G-19510** **EYE SCISSORS**  
angled to side, probe tips  
blunt tips  
overall length 115 mm



**G-31315 OPERATING SCISSORS**

straight blades  
sharp-blunt  
overall length 145 mm



**ENUCLEATION SCISSORS**

curved blades, blunt  
40 mm from tip to pivot  
overall length 135 mm

**G-19520**

strong curve

**G-19525**

medium curve

**G-19530**

light curve




UTILITY · ENUCLEATION · RING HANDLE




**20 GAUGE**

**VITREORETINAL SCISSORS**


20 gauge / 0.9 mm  
titanium handle  
gold-colored  
overall length 155 mm




45° angle  
horizontal  
**G-36502** 6 mm handle  
**G-36504** 8 mm handle




curved  
with membrane peeler  
**G-36506** 6 mm handle  
**G-36508** 8 mm handle




curved  
**G-36512** 6 mm handle  
**G-36514** 8 mm handle




*KROLL*  
angled, with membrane peeler  
**G-36516** 6 mm handle  
**G-36518** 8 mm handle




angled, blunt tips  
for MIS  
**G-36522** 6 mm handle  
**G-36524** 8 mm handle




*HEIMANN*  
angled  
**G-36526** 6 mm handle  
**G-36528** 8 mm handle




*HEIMANN*  
straight  
**G-36532** 6 mm handle  
**G-36534** 8 mm handle



*HEIMANN*  
straight, blunt tips  
**G-36536** 6 mm handle  
**G-36538** 8 mm handle




60° angled  
**G-36546** 6 mm handle  
**G-36548** 8 mm handle



curved  
horizontal  
**G-36552** 6 mm handle  
**G-36554** 8 mm handle



*NAITO*  
long blade 3.5 mm  
curved, tapered tip  
**G-36556** 6 mm handle  
**G-36558** 8 mm handle




*ROIDER*  
for RPE transplantation, vitrectomies  
and retinectomies  
**G-36572** 6 mm handle  
**G-36574** 8 mm handle


23 GAUGE

VITREORETINAL SCISSORS

23 gauge / 0.6 mm  
titanium handle  
gold-colored  
overall length 155 mm




curved



HEIMANN  
straight


**G-36566** 6 mm handle  
**G-36568** 8 mm handle

**G-36576** 6 mm handle  
**G-36578** 8 mm handle



HEIMANN  
for use with 23 gauge trocar system  
angled, sharp tips

**G-36540** 6 mm handle  
**G-36541** 8 mm handle



HEIMANN  
straight, blunt tips


**G-36582** 6 mm handle  
**G-36584** 8 mm handle

VITREORETINAL 23 & 25 GAUGE


25 GAUGE

VITREORETINAL SCISSORS

25 gauge / 0.5 mm  
titanium handle  
gold-colored



curved



HEIMANN  
straight

**G-36562** 6 mm handle  
**G-36564** 8 mm handle

**G-36586** 6 mm handle  
**G-36588** 8 mm handle


**ERGOLINE**

**ERGOLINE SCISSORS**


360° rotatable and interchangeable tips  
for Ergoline handles 6 and 8 mm diameter  
overall length 55 mm

**ERGOLINE HANDPIECE**

for Ergoline scissors  
gold-colored  
overall length 105 mm  
**G-36900** 6 mm titanium handpiece  
**G-36901** 8 mm titanium handpiece




**ERGOLINE SCISSORS**  
45° angled  
horizontal  
20 gauge / 0.9 mm  
**G-36902**




*KROLL*  
**ERGOLINE SCISSORS**  
angled, with membrane peeler  
20 gauge / 0.9 mm  
**G-36915**



*HEIMANN*  
**ERGOLINE SCISSORS**  
angled  
20 gauge / 0.9 mm  
**G-36925**



*HEIMANN*  
**ERGOLINE SCISSORS**  
straight  
**G-36930** 20 gauge / 0.9 mm  
**G-36975** 23 gauge / 0.6 mm  
**G-36985** 25 gauge / 0.5 mm



*NAITO*  
**ERGOLINE SCISSORS**  
long blade 3,5 mm, curved  
tip tapered to 0,75 mm  
20 gauge / 0.9 mm  
**G-36955** 20 gauge / 0.9 mm



**CLEANING ADAPTER**  
for rotatable and  
interchangeable tips  
**G-36604**

**09**

VITREORETINAL (INTERCHANGEABLE AND ROTATABLE)

SINGLE-USE  
INTERCHANGEABLE SCISSORS  
AVAILABLE ON P. 218



