






Quality Assurance



We at **Eagle Labs** guarantee that our products are manufactured with the highest quality standards. Eagle Labs complies with CFR 21 of The U.S. Food and Drug Administration. Our manufacturing facility is ISO 13485 and MDD 93/42/EEC Registered. Eagle Products carry the European Union CE Mark Registration.

Commitment

Our commitment at Eagle Labs remains the same as it has been for over 30 years. We continually strive to offer high-quality and cost-effective Single-Use Medical Devices to the Ophthalmic Community. We proudly distribute our products for surgery throughout the world and support those healthcare professionals who have dedicated their time and efforts to improve quality of sight and life for everyone.

Table of Contents

Anesthesia Administration

Atkinson Tip Peribulbar Needles	4
Atkinson Tip Retrobulbar Needles	4
Sharp Point Retrobulbar Needles	4
Tri-Port Sub-Tenon Anesthesia Cannulae	4

Anterior Chamber Maintenance

Anterior Chamber Maintainer	4
Anterior Chamber Maintainer Kit	4
Anterior Chamber Maintainer, Pediatric	4

Capsule Polishing

Aspirating J-Hook Capsule Polisher	5
Capsule Polisher (Kratz)	5
Micro Loop Curette Polisher	5
Olive Tip Capsule Polisher	5
Silicone "Squeegee" Capsule Polishers	5

Capsulotomy

Formed Irrigating Cystotomes	6
Straight Irrigating Cystotomes	6
Angled Irrigating Cystotomes	7
Capsulorhexis Suction Cystotome (Dahan)	7
Curved Irrigating Cystotomes	7

Corneal Abrasion

Stromal Puncture Needle	7
-----------------------------------	---

Cortex Removal

Aspirating J-Hook Cannulae	8
Aspirating J-Hook Capsule Polisher	8
Aspirating C-Hook Cannulae	8
Aspirating J-Hook "Squeegee" Cannulae	8
Bimanual I/A Cortex Removal System	8
Cortex Aspiration Cannulae	8
Simcoe Cortex Extraction Cannulae	8

Drapes & Fluid Collection

Small Plastic Drapes	9
Refractive Drapes	9
Medium Plastic Drapes	9
Half Body Drapes	10
Phaco Drainage Bags	10
Micro Surgical Accessory Kit	10
Fluid Drainage Wick	10
Ophthalmic Surgery Kit	10

Eye Protection

Bilateral Eye Shields, Clear	11
Corneal Light Shields	11
Eye Shields	11
Universal Eye Shields, Clear	11

FLACS

FLACS Cannula	11
-------------------------	----

Fluid Administration

Ocular Irrigation Tube	12
Ocular Irrigation Set	12
High Flow Filter	12

Table of Contents

Glaucoma

Anterior Chamber Cannulae	13
Glaucoma Revision Picks	13
Visco Canalostomy Cannulae	13

Hydrodissection

Cortical Cleaving Hydrodissection Cannulae	13, 14
Nucleus Hydrodissection Cannulae	13, 14
Chang Hydrodissection Cannulae	14

Hydrodelineation

Nucleus Hydrodelineation Cannulae	15
---	----

Irrigation, Air Injection & Viscoelastic

Anterior Chamber Washout Cannulae	15
Bishop Harmon Cannulae	15
Anterior Chamber Irrigation Cannulae	15, 16
FLOvisc Viscoelastic Cannulae	17
Stromal Hydration Cannulae	17
Infusion Cannulae	17

Irrigation & Aspiration

Double I and A Cannulae	18
Double I and A Simcoe Cannulae	18

Manipulators **GOLD™**

IOL Manipulators	30
----------------------------	----

Nucleus Removal

Hydro-Expression Cannulae	19
Irrigating Vectis	19
Viscoelastic Expression Cannulae	19

Oculoplastic

Lacrimal Cannulae	19
Lacrimal Intubation Sets	20

Refractive (LASIK)

LASIK Irrigating Cannulae	20-23
PVA LASIK Drainage Ring (Chayet)	23
LASIK Drainage Ring	23
LASIK Sponge	23
LASIK Shields	23
PVA LASIK Eye Spears	24
Pre-Expanded LASIK Spears	24
Refractive Drapes	24

Speculums

Wire Eye Speculums	25
------------------------------	----

Surgical Absorption Products

Cellulose Eye Sponges	25
Diamond Knife Cleaning Block	25
Eye Pledget	25
Instrument Wipes	25
Pre-Expanded LASIK Spears	25
PVA Eye Spears	25

Table of Contents

Surgical Knives **GOLD™**

Sub 2.0mm Phaco Slit Knives	31
Double Bevel Sideport Knives	31
Spoon Knives	31
MVR Knives	31
MVR Knives, Angled	32
Clear Corneal Knives	32
Trapezoid Knives	32
Stab Knives	33
Phaco Slit Knives	33
IOL Inserter Knives	33
Crescent Knives	33

Surgical Knives - Guarded **GOLD™**

Instructions For Use	34
Sub 2.0mm Phaco Slit Knives - Guarded	35
Clear Corneal Knives - Guarded	35
Trapezoid Knives - Guarded	35
Stab Knives - Guarded	36
Phaco Slit Knives - Guarded	36
IOL Inserter Knives - Guarded	36
Crescent Knives - Guarded	36

Surgical Blades

Microsurgical Blades, Depth Control	37
Needle Blades	38
Lamellar Blades	38
Scleral Blades	38
Miniature Edged Blades	38
Corneal Trephine Blades	39
Corneal Trephine Block	39
Ear Blades	39
Blade Holders	39

Syringes, Handles & Adapters

Spring Assisted Aspirating Syringes	26
Aspirating Syringes	26
Male/Male Luer Adapter	26
Irrigating Handle	26
Aspirating Handle	26

Vitreoretinal

Retinal Dissecting Cannulae	26
Retinal Micro Pick Cannulae	26
Hard Tapered Extrusion Cannulae	27
Silicone Oil Cannulae	27
Silicone Oil Cannulae, Manual Delivery	27
Infusion Cannulae	27
High Viscosity Infusion Cannulae	27
Subretinal Fluid Cannula	28
Subretinal Injection Cannulae	28
Subretinal Aspiration Cannulae	29
Perfluorocarbon Cannulae	29
Vitreoretinal Viscodissector	29

Ordering Information	Inside Back Cover
----------------------	-------------------

Product Index	40-43
---------------	-------

Worldwide Distribution	44
------------------------	----

Anesthesia Administration

Atkinson Tip Retrobulbar

10 per box

101-23	23GA x 1 1/2" (38mm)
101-25	25GA x 1 1/2" (38mm)
101-27	27GA x 1 1/2" (38mm)

Atkinson Tip Peribulbar

10 per box

102-23-32	23GA x 1 1/4" (32mm)
102-25-32	25GA x 1 1/4" (32mm)
102-27-32	27GA x 1 1/4" (32mm)
102-23-22	23GA x 7/8" (22mm)
102-25-22	25GA x 7/8" (22mm)
102-27-22	27GA x 7/8" (22mm)

Tri-Port Sub-Tenon Cannula

10 per box

2 - .4mm side ports, open-end port

100-19	19GA x 1" angled (25mm)
100-19C	19GA x 1" curved (25mm)
100-21	21GA x 1" angled (25mm)
100-21C	21GA x 1" curved (25mm)

Sharp Point Retrobulbar

10 per box

103-23	23GA x 1 1/2" (38mm)
103-25	25GA x 1 1/2" (38mm)

Anterior Chamber Maintenance

Anterior Chamber Maintainer Kit

3 per box

170-20T 20GA
(infusion cannula with self retaining serrated tip)
20GA x 3/16" (4.76mm)
cannula attached to 8" (20.3cm) flexible tubing
attached to an ocular irrigation set.

Anterior Chamber Maintainer

3 per box


170-20 20GA
(infusion cannula with self retaining serrated tip)
20GA x 3/16" (4.76mm) cannula
attached to 8" (20.3cm) flexible tubing

Pediatric A.C. Maintainer

3 per box

Used to maintain the chamber during pediatric retinal procedure (Mendez/Freeman)
171-20R 20GA
3.5mm, 90° angled, 20GA cannula
with 40° beveled tip connected to 11" (28cm)
flexible tube (.8mm I.D., 1.6mm O.D.)


Capsule Polishing



Silicone "Squeegee" Tip Capsule Polisher (Engels)

10 per box


40 degree curve with extended silicone sleeve for polishing
128-27 27GA x 1" (25mm)



Silicone "Squeegee" Tip Capsule Polisher (Terry)

10 per box


7mm bend to tip with flexible sleeve
125-23 23GA x 1" (25mm)
125-27 27GA x 1" (25mm)



Silicone "Squeegee" Tip Cortex Extractor and Polisher (Buedding)

10 per box


Irrigating cortex extractor with flexible sleeve for the 12 o'clock position
 straight 23mm long
107-27J 27GA



Capsule Polisher (Kratz)

10 per box


8mm bend, .30mm top port, rounded tip
125-21 21GA x 1" (25mm)
125-24 24GA x 1" (25mm)
 Roughened end permits scraping and aspiration of cortical debris



Capsule Polisher (Kratz)

10 per box


Roughened tip, angled 3mm
125-22 22GA x 1 1/2" (38mm)



Capsule Polisher (Lu)

10 per box



8mm bend, .30mm top port, rounded tip
124-23 23GA x 1" (25mm)
 Curved, roughened olive tip - open port



Micro Loop Curette Polisher

10 per box

125-25 25GA x 1 1/4" (32mm)
 Irrigating posterior capsule polisher
 curette loop is angled, posterior surface of loop is roughened, with open port











Olive Tip Capsule Polisher

10 per box

10mm bend with roughened tip - open port
126-23 23GA x 1" (25mm)
126-25 25GA x 1" (25mm)

Capsulotomy



Irrigating Cystotome

10 per box

13mm bend to tip (formed)

120-01 25GA x 3/4" (19mm)

120-02 27GA x 3/4" (19mm)

Irrigating Cystotome/Push-Pull

10 per box

13mm bend to tip (formed)

120-03 25GA x 3/4" (19mm)

Irrigating Cystotome

10 per box

Straight

120-21S 21GA x 1" (25mm)

120-22S 22GA x 1" (25mm)

120-23S 23GA x 1" (25mm)

120-25S 25GA x 5/8" (16mm)

120-27S 27GA x 1/2" (13mm)

120-30S 30GA x 3/4" (19mm)

Irrigating Cystotome

10 per box

15mm bend to tip (formed)

120-21F 21GA x 1" (25mm)

120-22F 22GA x 1" (25mm)

120-23F 23GA x 1" (25mm)

120-25F 25GA x 1" (25mm)

Irrigating Cystotome

10 per box

12mm bend to tip (formed)

120-22F-12 22GA x 3/4" (19mm)

120-23F-12 23GA x 3/4" (19mm)

120-25F-12 25GA x 5/8" (16mm)

120-27F-12 27GA x 1/2" (13mm)

10mm bend to tip (formed)

120-30F-10 30GA x 1/2" (13mm)

Irrigating Cystotome

10 per box

8mm bend to tip (formed)

120-25F-8 25GA x 5/8" (16mm)

120-27F-8 27GA x 1/2" (13mm)

Irrigating Cystotome/Push-Pull (Blumenthal)

10 per box

Straight

121-25S 25GA x 5/8" (16mm)

121-27S 27GA x 1/2" (13mm)

8mm bend to tip (formed)

121-25F-8 25GA x 5/8" (16mm)

121-27F-8 27GA x 1/2" (13mm)

12mm bend to tip (formed)

121-25F-12 25GA x 5/8" (16mm)

121-27F-12 27GA x 1/2" (13mm)

Capsulotomy

Irrigating Cystotome/Reverse Bend

10 per box

Straight

122-23S 23GA x 3/4" (19mm)

122-25S 25GA x 5/8" (16mm)

122-27S 27GA x 1/2" (13mm)

12mm bend to tip (formed)

122-23RF 23GA x 3/4" (19mm)

122-25RF 25GA x 5/8" (16mm)

122-27RF 27GA x 1/2" (13mm)

Irrigating Cystotome (Joyce)

10 per box

Curved for intercapsular/endocapsular techniques

220-25 25GA x 5/8" (16mm)

Irrigating Cystotome (Joyce)

Curved for intercapsular/endocapsular techniques

220-27 27GA x 1/2" (13mm)

Capsulorhexis Suction Cystotome (Dahan)

10 per box

Curved design to create suction
and hold capsular flap

221-25* 25GA x 1" (25mm)

*Used with extension tubing #169-30L-10
and syringe part #300-03 or #300-05

Irrigating Cystotome/Angled

10 per box

Angled 45°, 10mm bend to tip

123-25 25GA x 5/8" (16mm)

123-27 27GA x 1/2" (13mm)

Corneal Abrasion

Stromal Puncture Needle

10 per box

To achieve anterior stromal puncture of the epithelium

165-25 25GA x 1" (25mm)


angled 30°, 19mm bend to tip

sharp cutting tip


Cortex Removal





Irrigation





Aspiration





Aspiration




Male and female connectors are compatible with I/A and Phaco systems. Allows cortex removal through paracentesis or Phaco incision while maintaining the chamber with irrigation.




Hook dimensions 2.5mm long x 2.0mm wide




Hook dimensions 2.75mm long x 2.5mm wide




Hook dimensions 2.0mm long x 2.0mm wide



Hook dimensions 1.5mm long x 1.5mm wide



Hook dimensions 1.5mm long x 2.0mm wide



Bimanual I/A Cortex Removal System

10 per box

192-23

Irrigation 21GA x 1" (25mm)

Two .4mm ports with female adapter

Aspiration 23GA x 1" (25mm)

One .4mm port .5mm from tip with male connector

Cortex Aspiration Cannula

10 per box

191-21

21GA x 1" (25mm) .4mm port with male adapter

191-23

23GA x 1" (25mm) .4mm port with male adapter

191-25

25GA x 1" (25mm) .4mm port with male adapter

Aspirating J-Hook Cannula

10 per box

104-25

25GA x 1" (25mm)

Aspirating J-Hook Cannula/ Capsule Polisher (Thornton)

10 per box

Open end port, roughened on the end of hook for capsule polishing

107-25

25GA x 1" (25mm)

Aspirating J-Hook Cannula

10 per box

105-22L

22GA x 1" (25mm), Left

105-22R

22GA x 1" (25mm), Right

105-22S

22GA x 1" (25mm), Straight

105-25L

25GA x 1" (25mm), Left

105-25R

25GA x 1" (25mm), Right

105-25S

25GA x 1" (25mm), Straight

105-27L

27GA x 1" (25mm), Left

105-27R

27GA x 1" (25mm), Right

105-27S

27GA x 1" (25mm), Straight

Aspirating Cortex Extractor

10 per box

For aspirating cortex between the 10 and 2 o'clock positions

106-23S

23GA x 1" (25mm), Straight

106-25S

25GA x 1" (25mm), Straight

106-27L

27GA x 1" (25mm), Left

106-27R

27GA x 1" (25mm), Right

106-27S

27GA x 1" (25mm), Straight

106-30L

30GA x 1" (25mm), Left

106-30R

30GA x 1" (25mm), Right

106-30S

30GA x 1" (25mm), Straight

Aspirating C-Hook Extractor

10 per box

106-27SC

27GA

Simcoe Cortex Extractor Aspiration Cannula

10 per box

Curved with .3mm inside port


108-23

23GA x 1" (25mm)

108-25

25GA x 1" (25mm)


Drapes & Fluid Collection



Small Plastic Drape w/ Adhesive Aperture 10 per box

810-20


Overall size	16" x 16" (41cm x 41cm)
Adhesive aperture area	5" x 5" (13cm x 13cm)
Aperture circle	2.5" (6cm)



Small Plastic Drape w/ Adhesive Incise 10 per box

810-21


Overall size	16" x 16" (41cm x 41cm)
Adhesive aperture area	5" x 5" (13cm x 13cm)



Refractive Drape w/ Split Aperture 10 per box

810-22


Overall size	16" x 11" (41cm x 28cm)
Adhesive aperture area	3" (7.6cm)
Aperture oval width	1.5" (3.8cm)



Medium Plastic Drape w/ Adhesive Aperture 5 per box

810-30


Overall size	51" x 32" (129cm x 81cm)
Adhesive aperture area	6" x 5" (15cm x 13cm)
Aperture oval	2" x 3" (5cm x 8cm)



Medium Plastic Drape w/ Adhesive Incise 5 per box

810-60

Overall size	48" x 50" (122cm x 127cm)
Adhesive incise area	5" x 5" (13cm x 13cm)




Medium Plastic Drape w/ Adhesive Incise 5 per box

810-61

Overall size	48" x 51" (122cm x 130cm)
Adhesive incise area	4" x 5" (10cm x 13cm)


Drapes & Fluid Collection



Non-Woven Half Body Drape *5 per box* w/ Adhesive Incise & Fluid Collection Pouch

810-54


Overall Size 60" x 77" (152cm x 196cm)
Adhesive incise area 5" x 4" (13cm x 10cm)
with velcro cord holder



Non-Woven Half Body Drape *5 per box* w/ Adhesive Aperture & Fluid Collection Pouch

810-64


Overall Size 61" x 77" (155cm x 196cm)
Adhesive incise area 4" x 3" (10cm x 8cm)
Aperture oval 1" x 2" (3cm x 5cm)
with velcro cord holder



Phaco Drainage Bag - 400ml *20 per box*

700-03


3.875" x 8.25"; 400ml capacity with two drainage wicks and adhesive backing to adhere to surgical drape



Drainage Bag - 80ml *20 per box*

700-02

2" x 8.25"; 80ml capacity with one drainage wick and adhesive backing to adhere to surgical drape



Microsurgical Accessory Kit *20 per box*

700-06

Cost Saving Convenience Kit includes one each:
400ml Drainage Bag with Two Wicks, Corneal Light Shield,
Instrument Wipe



Fluid Drainage Wick *20 per box*

700-04

1/8" x 6.125" The biocompatible polymer is extremely hydrophilic and fast wicking, absorbing up to 25x its weight in fluids



Ophthalmic Surgery Kit *20 per box*

700-07

Cost Saving Convenience Kit includes one each:
Wick, Corneal Light Shield, Instrument Wipe

Eye Protection



Corneal Light Shield

20 per box

700-01

8mm diameter. Protects the retina from light. Contains no cellulose and leaves no lint or debris. Upon moistening, becomes soft and porous, conforming to shape of cornea. Keeps cornea moist. Easy to remove.



Universal Eye Shields, Clear

50 per box

400-15

Clear polycarbonate lenses, vented, non-sterile



Bilateral Eye Shield, Clear

25 per box

400-25

Clear polycarbonate lenses, vented, w/foam lining & adjustable strap, non-sterile




Eagle Eye Shield, White

25 per box

400-01

Eye shield, non-sterile, left or right eye




Eagle Eye Shield, White

10 per box

400-01S

Eye shield, sterile, left or right eye

FLACS



FLACS Cannula


10 per box

116-25

25GA x 1" (25mm)

10mm angled, spatula tip, used for dissection during femto cataract surgery

Fluid Administration




Ocular Irrigation Tube 20 per box

Extension tube with male/female luer connections provide fluid administration during ophthalmic procedures

169-30L-6 Irrigation tube - 6" (15.2cm)
(I.D. - .8mm, O.D. - 1.6mm)

169-30L-10 Irrigation tube - 10" (25.4cm)
(I.D. - .8mm, O.D. - 1.6mm)

169-30L-30 Irrigation tube - 30" (76.2cm)
(I.D. - 3mm, O.D. - 4mm)




Ocular Irrigation Set 3 per box


Features an adjustable roller clamp for control of fluid flow

Single
200-01 Tubing length 86"

Double
200-02 Tubing length 96"



Single w/ .20 micron filter
200-21 Tubing length 86"



Double w/ .20 micron filter
200-22 Tubing length 96"



High Flow Filter .20 Micron 50 per box

201-20
Male/female luer connections, used between syringe and cannula to filter out particulates and contaminants

Glaucoma

Visco Canalostomy Cannula

5 per box

177-33 5mm beveled tip, slightly angled
Micro-gauge cannula for viscoelastic injection during glaucoma surgery

Glaucoma Revision Pick

10 per box

60 Degree (Lederer)

177-27 27GA x 1/2" (13mm)
Angled 60°, 7mm bend to tip
Tip has a sharp side cutting edge (superior approach)

Glaucoma Revision Pick

10 per box

30 Degree (Lederer)

178-27 27GA x 1/2" (13mm)
Angled 30°, 10mm bend to tip
Tip has a sharp side cutting edge (temporal or nasal approach)

Anterior Chamber Cannula (Mendez)

10 per box

110-27TM 27GA x 10mm, round tip
Angled 45°, 3mm bend to tip to reform the anterior chamber during glaucoma procedures

Hydrodissection

Nucleus Hydrodissection Cannula

10 per box

1" (25mm) long with flattened tip

113-23 23GA 7mm bend to tip
113-25 25GA 7mm bend to tip
113-27 27GA 7mm bend to tip
113-30 30GA 7mm bend to tip
113-25-8 25GA 8mm bend to tip
113-23-10 23GA 10mm bend to tip
113-25-10 25GA 10mm bend to tip
113-27-10 27GA 10mm bend to tip
113-30-10 30GA 10mm bend to tip

Hydrodissection Cannula (Tuebingen)

10 per box

313-27 27GA x 1" (25mm), 75° curved, round tip
313-30 30GA x 1" (25mm), 75° curved, round tip

Nucleus Hydrodissection Cannula

10 per box

1" (25mm) long with flattened tip, vertical flow

213-25 25GA 7mm bend to tip
213-27 27GA 7mm bend to tip


Nucleus Hydrodissection Cannula (Nichamin)

10 per box

1" (25mm) long with flattened tip, formed

113-25F-12 25GA 12mm bend to tip
113-27F-12 27GA 12mm bend to tip
113-30F-12 30GA 12mm bend to tip

Hydrodissection




Chang Hydrodissection Cannula

10 per box

Angled 90°, 1.5mm flat beveled tip

113-01 27GA x 1" (25mm)




Nucleus Hydrodissector (Bobbitt)

10 per box

Double angled cannula allows insertion into anterior capsule for hydrodissection

116-27 27GA x 1" (25mm)




Cortical Cleaving Hydrodissector

10 per box

Angled 45°, 1" (25mm) long with round tip

216-30 30GA 12mm bend to tip




Cortical Cleaving Hydrodissector (Grabow)

10 per box

Angled 45°, 1" (25mm) long with round tip

217-27 27GA 11mm bend to tip




Cortical Cleaving Hydrodissector

10 per box

Angled 45°, 1" (25mm) long with flattened tip

218-27 27GA 11mm bend to tip




J-Hook Hydrodissector

10 per box

For hydrodissection between the 10 - 2 o'clock positions

105-27FL 27GA

Hook dimensions 2.25mm long X 1.75mm wide, flattened tip



Nucleus Hydrodissector (Johnson)


10 per box

115-25 25GA x 1" (25mm)


115-27 27GA x 1" (25mm)

115-30 30GA x 1" (25mm)

10mm bend with beveled tip for ease of insertion thru small incision.



horizontal flow



vertical flow

Nucleus Hydrodissector

10 per box


118-25 25GA x 1" (25mm) curved (Helsinki) horizontal

118-25S 25GA x 1" (25mm) curved (Sauter) vertical

118-27 27GA x 1" (25mm) curved (Helsinki) horizontal

118-27S 27GA x 1" (25mm) curved (Sauter) vertical

Hydrodelineation



Nucleus Hydrodelineation (Johnson)


10 per box

114-25 25GA x 1" (25mm)

114-27 27GA x 1" (25mm)

114-30 30GA x 1" (25mm)

With beveled tip for ease of insertion thru small incision




Nucleus Hydrodelineation

10 per box

117-23 23GA

117-25 25GA

Parallel bevel of cannula allows hydrodynamic penetration of the cataract, separating the lenticular lamellae of the inner nucleus, without plugging the cannula lumen




Nucleus Hydrodelineator/Rotator

10 per box

Formed 12mm from bend to tip. Tip is single bevel, angled 90°

119-25 25GA x 1" (25mm) curved

Irrigation, Air Injection & Viscoelastic




Anterior Chamber Washout Cannula

1" (25mm) long with flattened tip

111-16 16GA


111-19 19GA



Bishop Harmon Cannula

111-16BH 16GA x 1" (25mm), 6mm bend to tip

111-19BH 19GA x 1" (25mm), 8mm bend to tip



Anterior Chamber Irrigator, Straight

1" (25mm) long with round tip


112-20 20GA

112-23 23GA

112-25 25GA

112-27 27GA

112-30 30GA



Anterior Chamber Irrigator, Curved

112-25C 25GA x 1" (25mm) curved with round tip

Irrigation, Air Injection & Viscoelastic

Anterior Chamber Irrigator, Angled

10 per box

4mm bend to tip, angled 30°

110-21-4	21GA x 1" (25mm)
110-23-4	23GA x 1" (25mm)
110-25-4	25GA x 1" (25mm)
110-27-4	27GA x 1" (25mm)
110-30-4	30GA x 1" (25mm)

7mm bend to tip, angled 30°

110-20	20GA x 1" (25mm)
110-21	21GA x 1" (25mm)
110-22	22GA x 1" (25mm)
110-23	23GA x 1" (25mm)
110-25	25GA x 1" (25mm)
110-26	26GA x 1" (25mm)
110-27	27GA x 1" (25mm)
110-30	30GA x 1" (25mm)

9mm bend to tip, angled 30°

110-25-9	25GA x 1" (25mm)
110-27-9	27GA x 1" (25mm)
110-30-9	30GA x 1" (25mm)

11mm bend to tip, angled 30°







110-25-11	25GA x 1" (25mm)
110-27-11	27GA x 1" (25mm)
110-30-11	30GA x 1" (25mm)

8mm bend to tip, angled 25°

110-26-1/2	26GA x 1/2" (12mm)
-------------------	--------------------

12mm bend with beveled tip, angled 45°

110-30SP	26GA x 9/16" (14mm)
-----------------	---------------------



Irrigation, Air Injection & Viscoelastic

FLOvisc Viscoelastic Cannula (DeCamp)

10 per box

109-2025 25GA tip
20GA cannula, stepped down to a 25GA cannula,
10mm bend to tip, angled 30°, 20GA cannula
Insert extends into the hub to further reduce bubbling

109-2027 27GA tip
20GA cannula, stepped down to a 27GA cannula,
10mm bend to tip, angled 30°, 20GA cannula
Insert extends into the hub to further reduce bubbling



FLOvisc Viscoelastic Cannula (DeCamp)

10 per box

109-22 27GA tip
22GA cannula, stepped down to a 27GA cannula,
10mm bend to tip, angled 30°, 27GA cannula
extends into the hub to further reduce bubbling

109-27 27GA tip
20GA cannula, stepped down to a 27GA cannula,
10mm bend to tip, angled 30°

109-27-25 25GA tip
20GA cannula, stepped down to a 25GA cannula,
10mm bend to tip, angled 30°



FLOvisc Cannulae are specifically designed for easier expression of viscoelastics. the internal insert eliminates any air pocket to reduce bubbling.

Stromal Hydration Cannula

10 per box

165-27 27GA x 1" (25mm)
angled, 12mm bend closed end flattened tip, inside port
Designed to irrigate the outer stroma
for improved sealing of clear corneal incisions

Infusion Cannula (Murphy)


10 per box

170-20M 20GA x 1" (25mm)
angled 45°, 7mm bend to tip, 2 EA.
.4mm side ports, 45° beveled tip (up)



Infusion Cannula (Wright)

10 per box

170-20W 20GA x 1" (25mm)
angled 45°, 7mm bend to tip, 2 EA.
.4mm side ports, 45° beveled tip (down)





Irrigation & Aspiration



I and A Cannula (Simcoe)

6 per box



150-23 23GA x 23GA
Aspiration through 13" (33cm) flexible tubing with .4mm side port and irrigation through cannula



Reverse I and A Cannula (Simcoe)

6 per box



150-23R 23GA x 23GA
Aspiration through cannula with .4mm side port, irrigation through 13" (33cm) flexible tubing



Double I and A Cannula

10 per box

130-25 23GA x 25GA x 1 1/2" (38mm)
angled 35°, 8mm from bend to tip





Double I and A Cannula

10 per box

130-19 21GA x 19GA x 1 1/2" (38mm)
Irrigation/aspiration cannula


130-21 21GA x 21GA x 1 1/2" (38mm)
Irrigation/aspiration cannula



Double I and A Cannula

10 per box


130-23 23GA x 23GA x 1 1/2" (38mm)
with 4mm angled aspiration tip



Double I and A Cannula

10 per box

130-23L 23GA x 23GA x 1 1/2" (38mm)
Allows continuous irrigation/aspiration with iris retraction



Double I and A Cannula

10 per box

130-23J 23GA x 23GA x 1 1/2" (38mm)
With J - hook aspiration tip

Nucleus Removal

Viscoelastic Expression Cannula (Corydon) 10 per box

109-25 25GA x 1" (25mm)
Curved design allows placement under the nucleus for viscoelastic expression during ECCE.
Angled tip can be used to hook and remove nucleus.

Hydro-Expression Cannula (Corydon) 10 per box

119-27 27GA x 1" (25mm)
Curved design allows placement under the nucleus for expression with balanced salt solution during ECCE.
Angled tip can be used to hook and remove nucleus.

Irrigating Vectis 10 per box

140-25-1 25GA, 5mm wide with (1 ea.) .3mm open port
140-25-2 25GA, 5mm wide with (2 ea.) .3mm open port
140-25-3 25GA, 5mm wide with (3 ea.) .3mm open port
140-27-1 27GA, 3mm wide with (1 ea.) .3mm open port

Oculoplastic

Lacrimal Cannula 10 per box

Open end
262-26 26GA x 1/2" (12mm) straight

260-26 26GA x 1" (25mm) curved

160-21 21GA x 1 1/2" (38mm) curved
160-25 25GA x 1 1/2" (38mm) curved
160-26 26GA x 1 1/2" (38mm) curved


Closed end with side port

161-23 23GA x 1 1/2" (38mm) curved, .3mm port
161-25 25GA x 1 1/2" (38mm) curved, .3mm port
163-21 21GA x 1 1/2" (38mm) straight, .3mm port
163-25 25GA x 1 1/2" (38mm) straight, .3mm port

Lacrimal Cannula (Straight) 10 per box

162-23 19GA x 1" (25mm) with 1/2" (13mm) 23GA tip

Oculoplastic




Lacrimal Intubation Set

3 per box

180-23 Probe 23GA x 7" (18cm)
12" (30.5cm) flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)



180-27 Probe 27GA x 4 1/2" (11.5cm)
12" flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)



Lacrimal Intubation Set, (Crawford)

3 per box

180-23T Probe 23GA x 7" (18cm), olive tip
12" (30.5cm) flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)



180-27T Probe 27GA x 4 1/2" (11.5cm), olive tip
12" (30.5cm) flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)




Lacrimal Intubation Set, (Guibor)

3 per box

180-23C Probe 23GA x 7" (18cm), curved
12" (30.5cm) flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)



Refractive / Lasik



Linn LASIK Cannula

10 per box



619-25 25GA x 1" (25mm) angled
9mm bend open end flattened tip, 2 ports



LASIK Irrigating Cannula (Banaji)

10 per box



604-25 25GA x 1" (25mm)
Formed 12mm bend to tip, 6 ports, closed end



LASIK Irrigating Cannula (Guell)

10 per box

605-25 25GA x 1" (25mm)
Formed 9mm bend to tip, 6 ports, open end




LASIK Irrigating Cannula (Burrato)

10 per box

606-25 25GA x 1" (25mm)
Formed 9mm bend to tip, 2 ports, open end

Refractive / Lasik





LASIK Irrigating Cannula (Pettigrove)

10 per box

607-25

25GA x 1" (25mm)

Formed 12mm bend to tip, 6 ports, open end





LASIK Irrigating Cannula (Zuckerman)

10 per box

608-25

25GA x 1" (25mm)

Angled 9mm bend to tip, 6 ports, open end





Irrigating LASIK Spatula (Slade)

10 per box

613-26

26GA x 1" (25mm)

Formed 9mm bend to tip, spatula tip





LASIK Flap Irrigator (Vidaurre)

10 per box

600-06

Dual 25GA multi-port Irrigation tips


Angled 13mm bend to tip. 4 ports each tip



600-16

Dual 25GA multi-port Irrigation tips

Angled 9mm bend to tip. 2 ports each tip



LASIK Flap Irrigator/Squeegee (Magruder)

10 per box


602-23

23GA x 1" (25mm) angled, 7mm bend to tip

602-27

27GA x 1" (25mm) angled, 7mm bend to tip

Angled silicone tip allows for safe lifting of the flap, silicone sleeve acts as a "squeegee" to remove excess fluid and smooth down the flap to realign gutter



LASIK Flap Spatula/Irrigator (Hersh)

10 per box

610-23


23GA x 1" (25mm) angled, 9mm bend to tip

610-25

25GA x 1" (25mm) angled, 9mm bend to tip

610-27

27GA x 1" (25mm) angled, 9mm bend to tip



611-23

23GA x 1" (25mm) angled, 4mm bend to tip

611-25


25GA x 1" (25mm) angled, 4mm bend to tip

611-27

27GA x 1" (25mm) angled, 4mm bend to tip

Non invasive spatula tip allows for safe lifting of the flap, curved angle for smoothing the flap during alignment, open end for irrigation


Refractive / Lasik



LASIK Irrigating Cannula (Lindstrom) **10 per box**

603-25 25GA x 1" (25mm) angled 7mm bend to tip

Features three each .3mm ports used to irrigate the anterior and posterior sides of the flap, smooth angle curvature allows for safe flap alignment




LASIK Fountain Cannula **10 per box**

600-21 25GA x 1" (25mm), right

600-22 25GA x 1" (25mm), left


Both fountain cannulae are 25GA with a single port and a 10mm curved angle to conform to the corneal shape. The tip is a smooth closed end design slightly tapered to allow insertion under the flap.



Quad Ported LASIK Irrigating Cannula **10 per box**

600-09 23GA x 1 1/2" (38mm), 4ea. .3mm ports

This instrument features four ports to irrigate the surgical interface and flush any debris collected during surgery




LASIK Irrigating Cannula (Nichamin) **10 per box**

600-23 23GA x 1" (25mm) 12mm bend to tip

600-25 25GA x 1" (25mm) 12mm bend to tip

600-27 27GA x 1" (25mm) 12mm bend to tip

This instrument features a flattened tip with open-end irrigation, formed curvature used for elevating and smoothing the corneal flap




LASIK Irrigating Cannula (Nichamin) **10 per box**

601-23 23GA x 1" (25mm) 6mm bend to tip

601-25 25GA x 1" (25mm) 6mm bend to tip

601-27 27GA x 1" (25mm) 6mm bend to tip


This instrument features a flattened tip with open-end irrigation, formed curvature used for elevating and smoothing the corneal flap



LASIK Aspiration Spoon **10 per box**

600-07 25GA multi-port irrigation tip


This unique spoon-shaped instrument removes fluid and debris from beneath the flap




LASIK Irrigating Cannula **10 per box**

612-25 25GA X 1" (25mm)

Angled 45°, 12mm bend to tip




Refractive / Lasik



LASIK Irrigating Cannula (Soloman)

10 per box


614-25 25GA X 1" (25mm)
Angled / formed spatula tip with 3 ports



LASIK Irrigating Cannula (Donnenfeld)

10 per box

615-23 23GA X 1" (25mm)
Curved with a horizontal flattened tip



LASIK Irrigating Cannula

10 per box

617-27 27GA X 1" (25mm)
Formed, 10mm bend to tip, flattened tip, 4 ports

LASIK Drains & Shields from Eagle Laboratories are made of the highest quality **LINT-FREE** PVA material to consistently ensure a particulate-free environment during LASIK Surgery. Drains are placed onto the orbital area after the microkeratome is used, before the ap is lifted. Used to control stromal bed hydration.



LASIK Drainage Ring (Chayet)

20 per box

600-20 (12mm opening)



LASIK Drainage Ring

20 per box

600-30 Superior Hinge, (9.5mm opening)



LASIK Sponge

20 per box

600-40 (11mm opening)

LASIK shields are placed over the posterior side of the ap and hinge, to moisten and protect against unwanted laser ablation.



LASIK Shield

20 per box

600-08 Nasal Hinge, 7mm dia., (5mm edge)



620-08 Nasal Hinge, 8mm dia., (6mm edge)



620-09 Nasal Hinge, 9mm Dia., (7mm edge)



620-10 Superior Hinge, Dual Edge, (4mm & 8mm edges)

Refractive / Lasik



PVA Eye Spear

10 packs of 5

600-10

- **LINT-FREE** PVA material.
- Compressed spears are ideal for controlling and absorbing fluid in and around the orbital area during **LASIK**.




Pre-Expanded LASIK Spear

20 packs of 2

600-15

- Pre-expanded PVA eye spears offer a **LINT-FREE** absorption surface.
- Faster wicking than compressed spears.
- Pre-expanded spears are ideal for aligning and smoothing the corneal ap.




Small Plastic Drape w/ Adhesive Aperture

10 per box

810-20

Overall size	16" x 16" (41cm x 41cm)
Adhesive area	5" x 5" (13cm x 13cm)
Aperture, circle	2.5" (6cm)




Small Plastic Drape w/ Adhesive Incise

10 per box

810-21

Overall size	16" x 16" (41cm x 41cm)
Adhesive area	5" x 5" (13cm x 13cm)




Refractive Drape w/ Split Aperture

10 per box

810-22

Overall size	16" x 11" (40cm x 28cm)
Adhesive Aperture Width	3" (7.6cm)
Aperture oval width	1.5" (3.8cm)

Speculums



Wire Eye Speculum (Barraquer)

3 per box

940-10

10mm

Wire Eye Speculum (Barraquer)

3 per box

940-15

15mm

Wire Eye Speculum (Kratz / Barraquer)

3 per box

940-20

15mm


Wire Eye Speculum (Kratz / Barraquer)

3 per box

940-29

15mm (temporal)

Surgical Absorption Products



Instrument Wipe

20 per box

700-05

3" x 3"; medium density. Useful for cleaning delicate surgical instruments without leaving lint or debris.



Diamond Knife Cleaning Block

10 per box

700-15

Small pore construction absorbs fluids and removes debris from diamond blade when gentle pressure is applied.



Cellulose Eye Sponge

10 packs of 10

500-01

High volume fast wicking action for absorption of blood and fluids within the operating field.



PVA Eye Spear

10 packs of 5

600-10

- **LINT FREE** PVA material
- Ideal for corneal and refractive surgery



Pre-Expanded LASIK Spear

20 packs of 2

600-15

- Pre-expanded PVA eye spears offer a **LINT-FREE** absorption surface.
- Faster wicking than compressed spears.
- Pre-expanded spears are ideal for aligning and smoothing the corneal flap.



Eye Pledget

10 packs of 10

500-05

2mm x 6mm

LINT FREE PVA Sponge

Syringes, Handles & Adapters



Spring Assisted Aspirating Syringe (Pallin) 10 per box

- 300-03** 3ml aspirating syringe with Luer Lok
300-05 5ml aspirating syringe with Luer Lok



Aspirating Syringe 10 per box

- 300-03SP** 3ml aspirating syringe with Luer Lok
300-05SP 5ml aspirating syringe with Luer Lok



Male/Male Luer Adapter 50 per box

- 191-01**
 Male luer fittings on both ends, used for connecting Phaco and I/A tubing to instruments.



Irrigating Handle 10 per box


- 191-02** 3" (7.5cm)
 Luer fittings both ends



Aspirating Handle 10 per box


- 191-03** 3" (7.5cm)
 Luer fittings both ends, with aspiration port.

Vitreoretinal




Membrane Dissector (Kingham) 10 per box

- 169-21** 21GA x 1" (25mm)
 Sharp tip with cutting edge angled 45°.




Vitreoretinal Micro Pick (Glaser) 10 per box

- 169-25** 20GA x 1" (25mm) with 25GA x 3/16" (4.76mm) tip
 Sharp tip with cutting edge angled 50°.




Membrane Micro Pick (Eaton) 10 per box

- 169-30** 20GA x 1" (25mm) with 30GA x 3/16" (4.76mm) tip
 Sharp tip with cutting edge angled 50°.



Epiretinal Membrane Dissector (Wagner) 10 per box


- 168-21** 21GA x 1" (25mm)
 Sharp tip with cutting edge angled 45°.



Vitreoretinal Micro Pick 10 per box


- 168-25** 20GA x 1" (25mm) with 25GA x 3/16" (4.76mm) tip
 Sharp tip with cutting edge angled 50°.

Vitreoretinal



Hard Tapered Extrusion Cannula (Tim Mendez) 10 per box

176-25 20GA x 1" with 25GA x 1mm tip



Silicone Oil Cannula 3 per box

167-20TS 20GA x 8mm

167-19TS 19GA x 8mm

167-16TS 16GA x 8mm, beveled tip


With 6" (15cm) high pressure extension tube pre-attached to a 10 cc ring control syringe.



Silicone Oil Cannula (Wagner) 3 per box

167-16T 16GA x 8mm with 6" (15cm)

High pressure extension tube.



Silicone Oil Cannula (Peyman) 10 per box

167-16 16GA x 8mm

167-19 19GA x 8mm

For injection or aspiration of high viscosity fluids in the repair of retinal detachments.



Silicone Oil Removal Kit (T. Mendez) 3 per box

164-19 19GA

80" PVC Tubing with connectors and 19GA straight, Cannula fits both the Alcon and Storz cassettes.

* Allows the safe removal of silicone oil in 10 minutes or less.



Infusion Cannula 3 per box

175-25 20GA X 2.5mm, beveled tip

175-40 20GA X 4mm, beveled tip


175-60 20GA X 6mm, beveled tip

High Viscosity Infusion Cannula 3 per box

175-14 20GA X 4mm, beveled tip

175-16 20GA X 6mm, beveled tip

Vitreoretinal




Subretinal Fluid Cannula 6mm (Grizzard) 10 per box

20GA cannula w/ silicone sleeve extended 6mm beyond end, overall length with tubing 37mm.

166-25 straight

166-25C curved

Soft silicone tip ensures safe entry through a macular tear or hole for aspiration of subretinal fluids.





Subretinal Fluid Cannula 3mm (Grizzard) 10 per box

20GA cannula w/ silicone sleeve extended 3mm beyond end, overall length with tubing 34mm.

166-20 straight

166-21 curved


Soft silicone tip ensures safe entry through a macular tear or hole for aspiration of subretinal fluids.



Vitreoretinal Injection Cannula Straight (Vander) 10 per box

172-30 20GA x 1" (25mm)


With 30GA x 5mm round tip.



Vitreoretinal Injection Cannula Angled (Vander) 10 per box

172-30C 20GA x 1" (25mm)


With 30GA round tip.



Micro Vitreoretinal Cannula Angled (Glaser) 10 per box

171-32 20GA x 1" (25mm)


With 32GA tip, 4mm long, angled 65° 1mm from tip.



Subretinal Injection Cannula Straight (Glaser) 10 per box

173-32 20GA x 1" (25mm)

With 32GAx 3mm, 35° beveled tip.



Subretinal Injection Cannula Curved (Glaser) 10 per box

173-32C 20GA x 1" (25mm)



With 32GAx 3mm, 35° beveled tip.

Vitreoretinal




Subretinal Aspiration Cannula (Landers) 10 per box

174-27 20GA x 1" (25mm)
With 27GA, 6.5mm tip, angled
2.5mm AT 60° from 22.5° beveled tip.



Vitreoretinal Aspiration Cannula (Landers) 10 per box

174-28 20GA x 1" (25mm)
With 27GA x 5/32" (4mm) round tip.




Subretinal Aspiration Cannula (Lambert) 10 per box

173-34 20GA x 1" (25mm)
With 32GA x 5/16" (8mm) curved beveled tip.



Vitreoretinal Injection Cannula (Mendez) 10 per box



175-32 20GA x 1 1/4" (32mm)
With 32GA (7mm) curved round tip.



Perfluorocarbon Co-axial I/A Cannula 10 per box

Used for perfluorocarbon liquid exchange during retinal procedures.

149-19 25GA X 1 1/2" (38mm)
With open end irrigation port
19GA outer dia. with .4mm aspiration port.



Perfluorocarbon I/A Cannula 10 per box

Used for perfluorocarbon liquid exchange during retinal procedures.

149-23 23GA X 1 1/2" (38mm)
With open end irrigation port
25GA aspiration cannula.

149-25 25GA X 1 1/2" (38mm)
With open end irrigation port
25GA aspiration cannula.




90°



70°



45°



Vitreoretinal Viscodissector (Glaser) 3 per box

169-30-90 20GA x 1" (25mm) with 30GA x 3/16" (4.76mm) tip 90°

169-30-70 20GA x 1" (25mm) with 30GA x 3/16" (4.76mm) tip 70°

169-30-45 20GA x 1" (25mm) with 30GA x 3/16" (4.76mm) tip 45°

30GA cannula angled 90°, 70° OR 45° with 1mm flattened tip. Cannula is attached to a Luer lock high pressure extension tube (.03" I.D. x 8" long).

Manipulators



IOL Manipulator (Jaffe/Maltzman)

6 per box

127-01

Fits 0.2mm diameter holes or larger, angled 10mm



IOL Manipulator

6 per box

127-02

Fits 0.2mm diameter holes or larger, angled 10mm



IOL Manipulator (Sinskey)

6 per box

127-03

Fits 0.2mm diameter holes or larger, straight



IOL Manipulator (Sinskey)

6 per box

127-04

Fits 0.2mm diameter holes or larger, angled 10mm



Button Tip Manipulator

6 per box

127-05

0.5mm button tip, angled 10mm

May also be used for iris retraction



IOL Manipulator (Fenzl)

6 per box

127-06

Fits 0.2mm diameter holes or larger, angled 10mm



Surgical Knives



Sub 2.0mm Phaco Slit Knife

6 per box

For bi-manual Phaco and micro-incision cataract procedures.

- EG-1014** 1.4mm, angled, full handle, bevel up
- EG-1015** 1.5mm, angled, full handle, bevel up
- EG-1016** 1.6mm, angled, full handle, bevel up
- EG-1017** 1.7mm, angled, full handle, bevel up
- EG-1018** 1.8mm, angled, full handle, bevel up
- EG-1020** 2.0mm, angled, full handle, bevel up



Double Bevel Sideport Knife

6 per box

To create a self-sealing incision during bi-manual Phaco procedures.

- EG-5010** 1.0mm, angled, full handle, double bevel
- EG-5012** 1.2mm, angled, full handle, double bevel
- EG-5014** 1.4mm, angled, full handle, double bevel



Spoon Knife

6 per box

For scleral tunnel incision.

- EG-6821** 3.0 mm, angled, full handle, double bevel



MVR Knife

6 per box

For vitreoretinal infusion port or paracentesis incision.

- EG-5555** 19GA, full handle
- EG-5560** 20GA, full handle
- EG-5563** 23GA, full handle
- EG-5565** 25GA, full handle

Surgical Knives



MVR Knife, Angled 6 per box

For vitreoretinal infusion port or paracentesis incision

- EG-5660** 20GA, full handle, angled
- EG-5663** 23GA, full handle, angled
- EG-5665** 25GA, full handle, angled



Clear Corneal Knife 6 per box

- EG-1122** 2.2 mm, angled, full handle, double bevel
- EG-1124** 2.4 mm, angled, full handle, double bevel
- EG-1127** 2.75 mm, angled, full handle, double bevel
- EG-1130** 3.0 mm, angled, full handle, double bevel
- EG-1132** 3.2 mm, angled, full handle, double bevel



Clear Corneal Knife, Parallel Sides 6 per box

- EG-1327** 2.75 mm, angled, full handle, double bevel
- EG-1330** 3.0 mm, angled, full handle, double bevel
- EG-1332** 3.2 mm, angled, full handle, double bevel



Trapezoid Clear Corneal / Implant Knife 6 per box

- EG-2832** 2.85 mm - 3.2 mm angled, full handle, bevel up
- EG-3032** 3.0 mm - 3.2 mm angled, full handle, bevel up



Dual-purpose design provides the ability to accomplish clear corneal and implant incisions with one knife.

Surgical Knives



Stab Knife 6 per box

For side port or scleral incision

- EG-1500** 15°, full handle
- EG-2250** 22.5°, full handle
- EG-3000** 30°, full handle
- EG-4500** 45°, full handle



Phaco Slit Knife 6 per box

For clear corneal or scleral incision

- EG-5510** 2.5 mm, angled, full handle, bevel up
- EG-5568** 2.65 mm, angled, full handle, bevel up
- EG-4738** 2.75 mm, angled, full handle, bevel up
- EG-4740** 2.85 mm, angled, full handle, bevel up
- EG-4741** 2.85 mm, straight, full handle
- EG-5570** 3.0 mm, angled, full handle, bevel up
- EG-5571** 3.0 mm, straight, full handle
- EG-5572** 3.2 mm, angled, full handle, bevel up
- EG-5573** 3.2 mm, straight, full handle
- EG-5580** 3.5 mm, angled, full handle, bevel up



Short Stop IOL Inserter Knife 6 per box

For IOL enlargement

- EG-9336** 3.5 mm, angled, full handle, bevel up
- EG-9340** 4.0 mm, angled, full handle, bevel up
- EG-9350** 5.2 mm, angled, full handle, bevel up




Crescent Knife 6 per box

For Scleral tunnel incision

- EG-9402** 2.0 mm, angled, full handle, bevel up


Instructions for using Guarded - Surgical Knives



Rotate the protective sleeve while gripping the end of the knife handle to release the “locking pin”.



Slide the protective sleeve back while gripping the end of the knife handle.



Rotate the protective sleeve while gripping the end of the knife handle to secure the “locking pin” into open position.



Ready for use.

Reverse the above steps to secure the protective cover in the closed position prior to disposal into an approved sharps container.

Guarded - Surgical Knives



Sub 2.0mm Phaco Slit Knife

6 per box

For bi-manual Phaco and micro-incision cataract procedures.

- EP-1014** 1.4mm, angled, full handle, bevel up
- EP-1015** 1.5mm, angled, full handle, bevel up
- EP-1016** 1.6mm, angled, full handle, bevel up
- EP-1017** 1.7mm, angled, full handle, bevel up
- EP-1018** 1.8mm, angled, full handle, bevel up
- EP-1020** 2.0mm, angled, full handle, bevel up



Clear Corneal Knife

6 per box

- EP-1122** 2.2 mm, angled, full handle, double bevel
- EP-1124** 2.4 mm, angled, full handle, double bevel
- EP-1127** 2.75 mm, angled, full handle, double bevel
- EP-1130** 3.0 mm, angled, full handle, double bevel
- EP-1132** 3.2 mm, angled, full handle, double bevel



Clear Corneal Knife, Parallel Sides

6 per box

- EP-1327** 2.75 mm, angled, full handle, double bevel
- EP-1330** 3.0 mm, angled, full handle, double bevel
- EP-1332** 3.2 mm, angled, full handle, double bevel



Trapezoid Clear Corneal / Implant Knife

6 per box

- EP-2832** 2.85 mm - 3.2 mm angled, full handle, bevel up
- EP-3032** 3.0 mm - 3.2 mm angled, full handle, bevel up

Dual-purpose design provides the ability to accomplish clear corneal and implant incisions with one knife.

Guarded - Surgical Knives



Stab Knife

6 per box

For side port or scleral incision

- EP-1500** 15°, full handle
- EP-2250** 22.5°, full handle
- EP-3000** 30°, full handle
- EP-4500** 45°, full handle



Phaco Slit Knife

6 per box

For clear corneal or scleral incision

- EP-5510** 2.5 mm, angled, full handle, bevel up
- EP-5568** 2.65 mm, angled, full handle, bevel up
- EP-4738** 2.75 mm, angled, full handle, bevel up
- EP-4740** 2.85 mm, angled, full handle, bevel up
- EP-4741** 2.85 mm, straight, full handle
- EP-5570** 3.0 mm, angled, full handle, bevel up
- EP-5571** 3.0 mm, straight, full handle
- EP-5572** 3.2 mm, angled, full handle, bevel up
- EP-5573** 3.2 mm, straight, full handle
- EP-5580** 3.5 mm, angled, full handle, bevel up



Short Stop IOL Inserter Knife

6 per box

For IOL enlargement

- EP-9336** 3.5 mm, angled, full handle, bevel up
- EP-9340** 4.0 mm, angled, full handle, bevel up
- EP-9350** 5.2 mm, angled, full handle, bevel up




Crescent Knife

6 per box

For Scleral tunnel incision

- EP-9402** 2.0 mm, angled, full handle, bevel up

Microsurgical Blades



Microsurgical Blade **6 per box**

EL-7511 1.5mm, 15°, Orange



EL-7512 1.5mm, 15°, Orange



Microsurgical Blade **6 per box**

EL-7513 3.0mm, 15°, Blue



EL-7514 3.0mm, 15°, Blue




Microsurgical Blade **6 per box**

EL-7515 5.0mm, 15°, Green



EL-7516 5.0mm, 15°, Green



Microsurgical Blade **6 per box**

EL-7522 4.0mm, 22.5°, Grey



EL-7523 4.0mm, 22.5°, Grey



Microsurgical Blade **6 per box**

EL-7530 3.5mm, 30°, Black



EL-7531 3.5mm, 30°, Black



Microsurgical Blade **6 per box**

EL-7545 3.0mm, 45°, Tan



EL-7546 3.0mm, 45°, Tan

Surgical Blades



Needle Blade *6 per box*
EL-5910 3.0mm



EL-5920 5.0mm



EL-5930 7.0mm



EL-5940 1.5mm



Lamellar Blade *6 per box*
EL-6600 60°



Scleral Blade *6 per box*
EL-6823H Full handle



Scleral Blade *6 per box*
EL-5700



Miniature Edged Blade *20 per box*
EL-6200 Chisel tip



EL-6400 Round tip



EL-6500 Lance tip



EL-6700 Curved tip



Miniature Edged Blade *6 per box*
EL-6300 Pointed tip



EL-6900 Sharp all around

Surgical Blades/Holders



Corneal Trephine Blade 1 per box

707-00	7.00mm
707-25	7.25mm
707-50	7.50mm
707-75	7.75mm
708-00	8.00mm
708-25	8.25mm
708-50	8.50mm
708-75	8.75mm



Corneal Trephine Block 1 per box

701-01

Ear Blade 6 per box

EL-7100	Myringotomy
EL-7120	Myringotomy spear cut
EL-7121	Myringotomy lance cut
EL-7200	Spoon 60°

Blade Holder 1 per box

EL-4312	7.5cm universal handle
EL-4310	10cm universal handle

Product Index

Product Number	Description	Box Quantity	Page	Product Number	Description	Box Quantity	Page
100-19	Tri-Port Sub-Tenon Cannula	10	4	110-26	Anterior Chamber Irrigator	10	16
100-19 C	Tri-Port Sub-Tenon Cannula	10	4	110-26-1/2	Anterior Chamber Irrigator	10	16
100-21	Tri-Port Sub-Tenon Cannula	10	4	110-27	Anterior Chamber Irrigator	10	16
100-21C	Tri-Port Sub-Tenon Cannula	10	4	110-27-9	Anterior Chamber Irrigator	10	16
101-23	Retrobulbar Atkinson Tip	10	4	110-27TM	Anterior Chamber Cannula (Mendez)	10	13
101-25	Retrobulbar Atkinson Tip	10	4	110-30	Anterior Chamber Irrigator	10	16
101-27	Retrobulbar Atkinson Tip	10	4	110-25-4	Anterior Chamber Irrigator	10	16
102-23-22	Peribulbar Atkinson Tip	10	4	110-27-4	Anterior Chamber Irrigator	10	16
102-25-22	Peribulbar Atkinson Tip	10	4	110-30-4	Anterior Chamber Irrigator	10	16
102-27-22	Peribulbar Atkinson Tip	10	4	110-25-9	Anterior Chamber Irrigator	10	16
102-23-32	Peribulbar Atkinson Tip	10	4	110-25-11	Anterior Chamber Irrigator	10	16
102-25-32	Peribulbar Atkinson Tip	10	4	110-27-11	Anterior Chamber Irrigator	10	16
102-27-32	Peribulbar Atkinson Tip	10	4	110-30-9	Anterior Chamber Irrigator	10	16
103-23	Retrobulbar Sharp Point	10	4	10-30-11	Anterior Chamber Irrigator	10	16
103-25	Retrobulbar Sharp Point	10	4	110-30SP	Anterior Chamber Irrigator	10	16
104-25	Aspirating J-Hook Cannula	10	8	111-16	Washout Cannula	10	15
105-22L	Aspirating J-Hook Cannula	10	8	111-19	Washout Cannula	10	15
105-22R	Aspirating J-Hook Cannula	10	8	111-16BH	Bishop Harmon Cannula	10	15
105-22S	Aspirating J-Hook Cannula	10	8	111-19BH	Bishop Harmon Cannula	10	15
105-25L	Aspirating J-Hook Cannula	10	8	112-20	Anterior Chamber Irrigator	10	15
105-25R	Aspirating J-Hook Cannula	10	8	112-23	Anterior Chamber Irrigator	10	15
105-25S	Aspirating J-Hook Cannula	10	8	112-25	Anterior Chamber Irrigator	10	15
105-27L	Aspirating J-Hook Cannula	10	8	112-25C	Anterior Chamber Irrigator, Curved	10	15
105-27R	Aspirating J-Hook Cannula	10	8	112-27	Anterior Chamber Irrigator	10	15
105-27S	Aspirating J-Hook Cannula	10	8	112-30	Anterior Chamber Irrigator	10	15
105-27FL	J-Hook Hydrodissection Cannula	10	14	113-01	Chang Hydrodissection Cannula	10	14
106-23S	Aspirating Cortex Extractor	10	8	113-23	Nucleus Hydrodissection Cannula	10	13
106-25S	Aspirating Cortex Extractor	10	8	113-25	Nucleus Hydrodissection Cannula	10	13
106-27L	Aspirating Cortex Extractor	10	8	113-27	Nucleus Hydrodissection Cannula	10	13
106-27R	Aspirating Cortex Extractor	10	8	113-30	Nucleus Hydrodissection Cannula	10	13
106-27S	Aspirating Cortex Extractor	10	8	113-25-8	Nucleus Hydrodissection Cannula	10	13
106-27SC	Aspirating Cortex Extractor	10	8	113-23-10	Nucleus Hydrodissection Cannula	10	13
106-30L	Aspirating Cortex Extractor	10	8	113-25-10	Nucleus Hydrodissection Cannula	10	13
106-30R	Aspirating Cortex Extractor	10	8	113-27-10	Nucleus Hydrodissection Cannula	10	13
106-30S	Aspirating Cortex Extractor	10	8	113-30-10	Nucleus Hydrodissection Cannula	10	13
107-25	Aspirating J-Hook Cannula	10	8	113-25F-12	Nucleus Hydrodissection Cannula	10	13
107-27J	Squeegee Cortex Extractor	10	5	113-27F-12	Nucleus Hydrodissection Cannula	10	13
108-23	Simcoe Cortex Extractor	10	8	113-30F-12	Nucleus Hydrodissection Cannula	10	13
108-25	Simcoe Cortex Extractor	10	8	114-25	Nucleus Hydrodelamination Cannula	10	15
109-20-25	Tri-Step Viscoelastic Cannula	10	17	114-27	Nucleus Hydrodelamination Cannula	10	15
109-20-27	Tri-Step Viscoelastic Cannula	10	17	114-30	Nucleus Hydrodelamination Cannula	10	15
109-22	FLOvisc Viscoelastic Cannula	10	17	115-25	Nucleus Hydrodissection Cannula	10	14
109-25	Viscoelastic Expression Cannula	10	19	115-27	Nucleus Hydrodissection Cannula	10	14
109-27	FLOvisc Viscoelastic Cannula	10	17	115-30	Nucleus Hydrodissection Cannula	10	14
109-27-25	FLOvisc Viscoelastic Cannula	10	17	116-25	FLACS Cannula	10	11
110-20	Anterior Chamber Irrigator	10	16	116-27	Nucleus Hydrodissection Cannula	10	14
110-21	Anterior Chamber Irrigator	10	16	117-23	Nucleus Hydrodelamination Cannula	10	15
110-21-4	Anterior Chamber Irrigator	10	16	117-25	Nucleus Hydrodelamination Cannula	10	15
110-22	Anterior Chamber Irrigator	10	16	118-25	Nucleus Hydrodissection Cannula	10	14
110-23	Anterior Chamber Irrigator	10	16	118-25S	Nucleus Hydrodissection Cannula	10	14
110-23-4	Anterior Chamber Irrigator	10	16	118-27	Nucleus Hydrodissection Cannula	10	14
110-25	Anterior Chamber Irrigator	10	16	118-27S	Nucleus Hydrodissection Cannula	10	14

Product Index

Product Number	Description	Box Quantity	Page	Product Number	Description	Box Quantity	Page
119-25	Nucleus Hydrodelineator/Rotator	10	15	130-19	Double I/A Cannula	10	18
119-27	Hydro-Expression Cannula	10	19	130-21	Double I/A Cannula	10	18
120-01	Irrigating Cystotome	10	6	130-23	Double I/A Cannula	10	18
120-02	Irrigating Cystotome	10	6	130-23J	Double I/A Cannula	10	18
120-03	Irrigating Cystotome	10	6	130-23L	Double I/A Cannula	10	18
120-21S	Irrigating Cystotome	10	6	130-25	Double I/A Cannula	10	18
120-22S	Irrigating Cystotome	10	6	140-25-1	Irrigating Vectis Loop	10	19
120-23S	Irrigating Cystotome	10	6	140-25-2	Irrigating Vectis Loop	10	19
120-25S	Irrigating Cystotome	10	6	140-25-3	Irrigating Vectis Loop	10	19
120-27S	Irrigating Cystotome	10	6	140-27-1	Irrigating Vectis Loop	10	19
120-30S	Irrigating Cystotome	10	6	149-19	Perfluorocarbon Co-axial I/A Cannula	10	29
120-21F	Irrigating Cystotome	10	6	149-23	Perfluorocarbon I/A Cannula	10	29
120-22F	Irrigating Cystotome	10	6	149-25	Perfluorocarbon I/A Cannula	10	29
120-23F	Irrigating Cystotome	10	6	150-23	Simcoe I/A Cannula	6	18
120-25F	Irrigating Cystotome	10	6	150-23R	Simcoe I/A Cannula Reverse	6	18
120-25F-8	Irrigating Cystotome	10	6	160-21	Lacrimonal Cannula	10	19
120-27F-8	Irrigating Cystotome	10	6	160-25	Lacrimonal Cannula	10	19
120-22F-12	Irrigating Cystotome	10	6	160-26	Lacrimonal Cannula	10	19
120-23F-12	Irrigating Cystotome	10	6	161-23	Lacrimonal Cannula	10	19
120-25F-12	Irrigating Cystotome	10	6	161-25	Lacrimonal Cannula	10	19
120-27F-12	Irrigating Cystotome	10	6	162-23	Lacrimonal Cannula	10	19
120-30F-10	Irrigating Cystotome	10	6	163-21	Lacrimonal Cannula	10	19
121-25S	Irrigating Cystotome/Push Pull	10	6	163-25	Lacrimonal Cannula	10	19
121-27S	Irrigating Cystotome/Push Pull	10	6	164-19	Silicone Oil Removal Kit	3	27
121-25F-8	Irrigating Cystotome/Push Pull	10	6	165-25	Stromal Puncture Needle	10	7
121-27F-8	Irrigating Cystotome/Push Pull	10	6	165-27	Stromal Hydration Cannula	10	17
121-25F-12	Irrigating Cystotome/Push Pull	10	6	166-20	Subretinal Fluid Cannula, 3mm	10	28
121-27F-12	Irrigating Cystotome/Push Pull	10	6	166-21	Subretinal Fluid Cannula, 3mm	10	28
122-23S	Irrigating Cystotome/Reverse Bend	10	7	166-25	Subretinal Fluid Cannula, 6mm	10	28
122-25S	Irrigating Cystotome/Reverse Bend	10	7	166-25C	Subretinal Fluid Cannula, 6mm	10	28
122-27S	Irrigating Cystotome/Reverse Bend	10	7	167-16	Silicone Oil Cannula	10	27
122-23RF	Irrigating Cystotome/Reverse Bend	10	7	167-16T	Silicone Oil Cannula	3	27
122-25RF	Irrigating Cystotome/Reverse Bend	10	7	167-16TS	Silicone Oil Cannula	3	27
122-27RF	Irrigating Cystotome/Reverse Bend	10	7	167-19	Silicone Oil Cannula	10	27
123-25	Irrigating Cystotome/Angled	10	7	167-19TS	Silicone Oil Cannula	3	27
123-27	Irrigating Cystotome/Angled	10	7	167-20TS	Silicone Oil Cannula	3	27
124-23	Capsule Polisher (Lu)	10	6	168-21	Membrane Dissector	10	27
125-21	Capsule Polisher	10	5	168-25	Vitreoretinal Micro Pick	10	26
125-22	Capsule Polisher	10	5	169-21	Membrane Dissector	10	26
125-23	Squeegee Capsule Polisher	10	5	169-25	Vitreoretinal Micro Pick	10	26
125-24	Capsule Polisher	10	5	169-30	Membrane Micro Pick	10	26
125-25	Micro Loop Polisher	10	5	169-30-90	Vitreoretinal Viscodissector	3	29
125-27	Squeegee Capsule Polisher	10	5	169-30-70	Vitreoretinal Viscodissector	3	29
126-23	Olive Tip Polisher	10	5	169-30-45	Vitreoretinal Viscodissector	3	29
126-25	Olive Tip Polisher	10	5	169-30L-6	Fluid Tubing Extension Set	20	12
127-01	IOL Manipulator	6	30	169-30L-10	Fluid Tubing Extension Set	20	12
127-02	IOL Manipulator	6	30	169-30L-30	Fluid Tubing Extension Set	20	12
127-03	IOL Manipulator	6	30	170-20	Anterior Chamber Maintainer	3	4
127-04	IOL Manipulator	6	30	170-20M	Infusion Cannula	10	17
127-05	Button Tip Manipulator	6	30	170-20T	Anterior Chamber Maintainer Kit	3	4
127-06	IOL Manipulator	6	30	170-20W	Infusion Cannula	10	17
128-27	Squeegee Capsule Polisher	10	5	171-20R	Anterior Chamber Maintainer	3	4

Product Index

Product Number	Description	Box Quantity	Page	Product Number	Description	Box Quantity	Page
171-32	Micro Vitreoretinal Cannula	10	28	3400-01S	Eye Shield	10	11
172-30	Vitreoretinal Injection Cannula	10	28	400-15	Universal Eye Shield, Clear (non-sterile)	50	11
172-30C	Vitreoretinal Injection Cannula	10	28	400-25	Bilateral Eye Shields, Clear	25	11
173-32	Subretinal Injection Cannula	10	28	500-01	Eye Sponges	10 packs of 10	25
173-32C	Subretinal Injection Cannula	10	28	500-05	Eye Pledget	10 packs of 10	25
173-34	Subretinal Injection Cannula (Lambert)	10	29	600-06	LASIK Flap Irrigator	10	20
174-27	Subretinal Injection Cannula	10	29	600-07	LASIK Aspiration Cannula	10	22
174-28	Vitreoretinal Aspiration Cannula (Landers)	10	29	600-08	LASIK Shield	20	23
175-14	High Viscosity Infusion Cannula	3	27	600-09	LASIK Irrigation Cannula	10	22
175-16	High Viscosity Infusion Cannula	3	27	600-10	PVA Eye Spears	10 packs of 5	24,25
175-25	Infusion Cannula 2.5mm	3	27	600-15	PVA Pre-Expanded LASIK Spears	20 packs of 2	24,25
175-32	Vitreoretinal Injection Cannula	10	29	600-16	LASIK Flap Irrigator	10	21
175-40	Infusion Cannula 4mm	3	27	600-20	LASIK Drainage Ring, (Chayet)	20	23
175-60	Infusion Cannula 6mm	3	27	600-21	LASIK Fountain Cannula, Right	10	22
176-25	Extrusion Cannula	10	27	600-22	LASIK Fountain Cannula, Left	10	22
177-27	Glaucoma Revision Pick	10	13	600-23	LASIK Irrigation Cannula	10	22
178-27	Glaucoma Revision Pick	10	13	600-25	LASIK Irrigation Cannula	10	22
177-33	Visco Canalostomy Cannula	5	13	600-27	LASIK Irrigation Cannula	10	22
180-23	Lacrimal Intubation Set	3	20	600-30	LASIK Drainage Ring (9.5mm opening)	20	23
180-23C	Lacrimal Intubation Set	3	20	600-40	LASIK Sponge	20	23
180-23T	Lacrimal Intubation Set	3	20	601-23	LASIK Irrigation Cannula	10	22
180-27	Lacrimal Intubation Set	3	20	601-25	LASIK Irrigation Cannula	10	22
180-27T	Lacrimal Intubation Set	3	20	601-27	LASIK Irrigation Cannula	10	22
191-01	Male/Male Luer Adapter	50	26	602-23	Silicone Tip LASIK Irrigator	10	21
191-02	Irrigating Handle	10	26	602-27	Silicone Tip LASIK Irrigator	10	21
191-03	Aspirating Handle	10	26	603-25	LASIK Irrigation Cannula	10	22
191-21	Cortex Aspiration Cannula	10	8	604-25	LASIK Irrigating Cannula	10	20
191-23	Cortex Aspiration Cannula	10	8	605-25	LASIK Irrigating Cannula	10	20
191-25	Cortex Aspiration Cannula	10	8	606-25	LASIK Irrigating Cannula	10	20
192-23	Bimanual I/A Cortex Removal System	10	8	607-25	LASIK Irrigating Cannula	10	20
200-01	Fluid Administration Sets Single	3	12	608-25	LASIK Irrigating Cannula	10	20
200-02	Fluid Administration Sets Double	3	12	610-23	LASIK Flap Spatula/Irrigator	10	21
200-21	Fluid Administration Sets Single Port w/Filter	3	12	610-25	LASIK Flap Spatula/Irrigator	10	21
200-22	Fluid Administration Sets Double Port w/Filter	3	12	610-27	LASIK Flap Spatula/Irrigator	10	21
201-20	High Flow Filter .20 Micron	50	12	611-23	LASIK Flap Spatula/Irrigator	10	21
213-25	Hydrodissection Cannula	10	13	611-25	LASIK Flap Spatula/Irrigator	10	21
213-27	Hydrodissection Cannula	10	13	611-27	LASIK Flap Spatula/Irrigator	10	21
216-30	Hydrodissection Cannula	10	14	612-25	LASIK Irrigating Cannula	10	22
217-27	Hydrodissection Cannula	10	14	613-26	Irrigating LASIK Spatula	10	20
218-27	Hydrodissection Cannula	10	14	614-25	LASIK Irrigating Cannula	10	23
220-25	Irrigating Cystotome	10	7	615-23	LASIK Irrigating Cannula	10	23
220-27	Irrigating Cystotome	10	7	617-27	LASIK Irrigating Cannula	10	23
221-25	Suction Cystotome	10	7	619-25	Linn LASIK Cannula	10	20
260-26	Lacrimal Cannula	10	19	620-08	LASIK Shield, 8mm Dia., (6mm edge)	20	23
262-26	Lacrimal Cannula	10	19	620-09	LASIK Shield, 9mm Dia., (7mm edge)	20	23
300-03	Spring Assisted Aspirating Syringe, 3ml	10	26	620-10	LASIK Shield, Dual Edge 4mm & 8mm	20	23
300-05	Spring Assisted Aspirating Syringe, 5ml	10	26	700-01	Corneal Light Shield	20	11
300-03SP	Aspirating Syringe, 3ml	10	26	700-02	Phaco Drainage Bag 80ml	20	10
300-05SP	Aspirating Syringe, 5ml	10	26	700-03	Phaco Drainage Bag 400ml	20	10
313-27	Hydrodissection Cannula (Tuebingen)	10	13	700-04	Fluid Drainage Wick	20	10
313-30	Hydrodissection Cannula (Tuebingen)	10	13	700-05	Instrument Wipe	20	25
400-01	Eye Shield, Non-Sterile	25	11	700-06	Microsurgical Accessory Kit	20	10

Product Index

Product Number	Description	Box Quantity	Page
700-07	Ophthalmic Surgery Kit	20	10
700-15	Diamond Knife Cleaning Block	10	25
701-01	Corneal Teflon Block	1	39
707-00	Corneal Trephine Blade	1	39
707-25	Corneal Trephine Blade	1	39
707-50	Corneal Trephine Blade	1	39
707-75	Corneal Trephine Blade	1	39
708-00	Corneal Trephine Blade	1	39
708-25	Corneal Trephine Blade	1	39
708-50	Corneal Trephine Blade	1	39
708-75	Corneal Trephine Blade	1	39
810-20	Small Plastic Drape w/Adhesive Aperture	10	9,24
810-21	Small Plastic Drape w/Adhesive Incise	10	9,24
810-22	Refractive Drape w/ Split Aperture	10	9,24
810-30	Medium Plastic Drape w/Aperture	5	9
810-54	Half Body Drape w/Incise	5	10
810-60	Medium Plastic Drape w/Incise	5	9
810-61	Medium Plastic Drape w/Incise	5	9
810-64	Half Body Drape w/Aperture	5	10
940-10	Wire Eye Speculum	3	25
940-15	Wire Eye Speculum	3	25
940-20	Wire Eye Speculum	3	25
940-29	Wire Eye Speculum	3	25
EG-1014	Sub 2.0mm Phaco Slit Knife	6	31
EG-1015	Sub 2.0mm Phaco Slit Knife	6	31
EG-1016	Sub 2.0mm Phaco Slit Knife	6	31
EG-1017	Sub 2.0mm Phaco Slit Knife	6	31
EG-1018	Sub 2.0mm Phaco Slit Knife	6	31
EG-1020	Sub 2.0mm Phaco Slit Knife	6	31
EG-1122	Clear Corneal Knife	6	32
EG-1124	Clear Corneal Knife	6	32
EG-1127	Clear Corneal Knife	6	32
EG-1130	Clear Corneal Knife	6	32
EG-1132	Clear Corneal Knife	6	32
EG-1327	Clear Corneal Knife, Parallel Sides	6	32
EG-1330	Clear Corneal Knife, Parallel Sides	6	32
EG-1332	Clear Corneal Knife, Parallel Sides	6	32
EG-1500	Stab Knife	6	33
EG-2250	Stab Knife	6	33
EG-2832	Trapezoid Clear Corneal/ Implant Knife	6	32
EG-3000	Stab Knife	6	33
EG-3032	Trapezoid Clear Corneal/ Implant Knife	6	32
EG-4500	Stab Knife	6	33
EG-4738	Phaco Slit Knife	6	33
EG-4740	Phaco Slit Knife	6	33
EG-4741	Phaco Slit Knife	6	33
EG-5010	Double Bevel Sideport Knife	6	31
EG-5012	Double Bevel Sideport Knife	6	31
EG-5014	Double Bevel Sideport Knife	6	31
EG-5510	Phaco Slit Knife	6	33
EG-5555	MVR Knife	6	31
EG-5560	MVR Knife	6	31

Product Number	Description	Box Quantity	Page
EG-5563	MVR Knife	6	31
EG-5565	MVR Knife	6	31
EG-5568	Phaco Slit Knife	6	33
EG-5570	Phaco Slit Knife	6	33
EG-5571	Phaco Slit Knife	6	33
EG-5572	Phaco Slit Knife	6	33
EG-5573	Phaco Slit Knife	6	33
EG-5580	Phaco Slit Knife	6	33
EG-5660	MVR Knife, Angled	6	32
EG-5663	MVR Knife, Angled	6	32
EG-5665	MVR Knife, Angled	6	32
EG-6821	Spoon Knife	6	31
EG-9336	Short Stop IOL Inserter Knife	6	33
EG-9340	Short Stop IOL Inserter Knife	6	33
EG-9350	Short Stop IOL Inserter Knife	6	33
EG-9402	Crescent Knife	6	33
EL-4310	10cm Universal Handle, Non-Sterile	1	39
EL-4312	7.5cm Universal Handle, Non-Sterile	1	39
EL-5700	Scleral Blade, CM	6	38
EL-5910	3.0mm Needle Blade, CM	6	38
EL-5920	5.0mm Needle Blade, CM	6	38
EL-5930	7.0mm Needle Blade, CM	6	38
EL-5940	1.5mm Needle Blade, CM	6	38
EL-6200	Chisel Tip Blade, CM	20	38
EL-6300	Pointed Tip Blade, CM	6	38
EL-6400	Round Tip Blade, CM	20	38
EL-6500	Lance Tip Blade, CM	20	38
EL-6600	60 Degree Lamellar Blade, CM	6	38
EL-6700	Curved Tip Blade, CM	20	38
EL-6823H	Scleral Blade, FH	6	38
EL-6900	Sharp Tip All Around, CM	6	38
EL-7100	Myringotomy Blade, CM	6	39
EL-7120	Myringotomy Blade, CM	6	39
EL-7121	Myringotomy Blade, CM	6	39
EL-7200	Myringotomy Blade, CM	6	39
EL-7511	1.5mm x 15 Degree Blade, CM	6	37
EL-7512	1.5mm x 15 Degree Blade, FH	6	37
EL-7513	3.0mm x 15 Degree Blade, CM	6	37
EL-7514	3.0mm x 15 Degree Blade, FH	6	37
EL-7515	5.0mm x 15 Degree Blade, CM	6	37
EL-7516	5.0mm x 15 Degree Blade, FH	6	37
EL-7522	4.0mm x 22.5 Degree Blade, CM	6	37
EL-7523	4.0mm x 22.5 Degree Blade, FH	6	37
EL-7530	3.5mm x 30 Degree Blade, CM	6	37
EL-7531	3.5mm x 30 Degree Blade, FH	6	37
EL-7545	3.0mm x 45 Degree Blade, CM	6	37
EL-7546	3.0mm x 45 Degree Blade, FH	6	37

Worldwide Distribution

Argentina

Rosinov Surgical
Instruments (GSJ S.A.)
www.moebiusqi.com.ar
info@rosinov.com
Phone: +54 11 48553668

Australia

Southern Cross
Hospital Supplies
mktg@schs.net.au
Phone: 02-9890-5088

Austria/

Eastern Europe

Askin & Co.
www.askin.co.at
sales@askin.co.at
Phone: +43/1/979 88 44

Belgium

DeCo-Vision
www.deco-vision.be
info@deco-vision.be
Phone: +32 (0) 9 245 47 72

Chile

Oftomed S.A.
www.oftomed.cl
mvillablanca@oftomed.cl
Phone: +56 2 571 39 90

Costa Rica

Insumed Inc.
www.insumedcr.com
jchacon@insumedcr.com
Phone: 506-2296-7422

Denmark

Lindcon Optical Group A/S
www.lindcon.dk
lindcon@lindcon.dk
Phone: 45 86 99 22 77

Egypt

Medi-Tec Egypt Corp
info@meditecegypt.com
Phone: 202 35717844

France

Carl Zeiss Meditec SAS
www.meditec.zeiss.com/iol
cataract-sales.meditec.fr@zeiss.com
+5-46-44-85-50

Germany

ascafamed GmbH
www.ascafamed.de
info@ascafamed.de
Phone: +49 (0) 6182. 92 92-10

Greece

Pharmex S.A.
www.pharmex.gr
info@pharmex.gr
Phone: +30 210 5199200

Indonesia

Pt. Universal Star Multilink
(USM)
sales.usm@gmail.com
Phone: (021) 2900 1517-19

Italy

Alfa Intes Casoria
www.alfaintes.it
info@alfaintes.it
Phone: +39 081 19527055

Japan

Nihon Surge
nihon-surge@nihon-surge.com
customer-support@nihon-surge.com
nihon-surge.com
Phone: +81 66 190 7865

Latvia

Latvijas-Amerikas
Acu Centrs
www.laac.lv
laac@laac.lv
Phone: +371 67272257

Mexico

Oftalmica Internacional
www.oftalmica.com

Tijuana Office

marcela@oftalmica.com
Phone: 01-800-026-6065

Mexico City Office

carmen@oftalmica.com
Phone: 01-800-015-2800

Norway

Medistim Norge AS
www.medistim.no
Orders:
bestillinger@medistim.com
Cust. Service:
norge@medistim.com
Phone: (+47) 23 03 52 50

Perú

Droguería ASG Medical
www.asgmedical.com.pe
ventas@asgmedical.com.pe
Phone: (511) 467-6767
(511) 467-6769

Poland

Consultronix S.A.
www.consultronix.pl
biuro@cxsa.pl
Phone: +48 12 290 22 22

Saudi Arabia

Arab Supply Est.
www.arabsupply.com.sa
ase@arabsupply.com.sa
Phone: 1-279-1526

South Africa

National Medical Supplies
www.nationalmedical.co.za
info@nationalmedical.co.za
Phone: 012-3865698

South Korea

Sung Woo Medical Co., Ltd.
sungwoomedi@daum.net
Phone: 2-783-4052

Spain

Bloss Group
www.blossgroup.com
bloss@blossgroup.com
Phone: +(34) 93 303 65 65

Sweden

Argus Medicinteknik AB
www.argusmedicinteknik.se
info@argusmedicinteknik.se
Phone: 0480 491280

Switzerland

medilas ag
www.medilas.ch
info@medilas.ch
Phone: +41 44 747 40 00

U.K. / Ireland

Dixey Instruments
www.dixeyinstruments.com
info@dixeyinstruments.com
Phone: +44 (0)1604 882480

United Kingdom

M.I.S.S. Ophthalmics Ltd.
www.miss-opth.com
missophthalmics@miss-opth.com
Phone: +44 (0)1536 741139

Uruguay

CLINIMED
www.clinimed.com.uy
clinimed@adinet.com.uy
infoclinimed@gmail.com
Phone: +598 23 36 05 67

Rosinov Surgical
Instruments (GSJ S.A.)
www.moebuisqi.com.ar
info@rosinov.com
Phone 598 2914 8875



How to Order

Toll Free Phone: 800 782-6534
Phone: 909 481-0011
Fax: 909 481-4481

Please contact us via phone, fax, email or mail and you will receive prompt assistance from a qualified customer service representative.

Order online: www.eaglelabs.com
by e-mail: sales@eaglelabs.com

or by US Mail: Eagle Labs
10201-A Trademark Street
Rancho Cucamonga, CA 91730-5850 USA

Shipping within U.S. – FOB Rancho Cucamonga, California. Shipping and insurance will be added to the order invoice.
Payment Terms – Net 30 days

Returned product must be accompanied by an RGA (Returned Goods Authorization) Number. Only merchandise credit will be issued unless a shipping error occurs from the manufacturer (Eagle Labs). Returns are subject to a 15% re-stocking fee. Product returns will be accepted within 6 months from ship date.

All Eagle Products are “Made in the USA” and unconditionally guaranteed against manufacturing defects.

Eagle products are single use and provided sterile unless otherwise noted.



Eagle Labs
10201-A Trademark Street
Rancho Cucamonga, CA 91730-5850 USA

Toll Free Ph 800 782-6534
Ph 909 481-0011
Fax 909 481-4481
goeagle@eaglelabs.com
www.eaglelabs.com

