

Knives

Tooke

Corneal Knife
straight blade
with curved cutting edge

K2-3650



Corneal Dissector

For intrastromal dissection.

K2-3660 angled, straight blade

K2-3661 angled, curved blade



K2-3660

K2-3661



Gill

Corneal Knife
curved cutting edge

K2-3700 straight

K2-3710 light curve



Fukasaku

"Hockey" Knife
sharp 3mm front edge
semi-sharp 6mm bottom edge

K2-3735



Spuds • Burrs

Foreign Body Spud

K2-4000 Golf Club, angled, lightly rounded tip

K2-4050 Francis, angled, flag shaped tip

K2-4100 Ellis, lightly curved and rounded tip



Spud & Gouge

in reversible handle

For removing foreign bodies.

K2-4600



Wills Eye Hospital

Foreign Body Spud & Burr
delicate curved spud on one
end and small burr for rust
ring removal on opposite end

K2-4800



Burrs

Algerbrush II

K2-4900 Corneal Rust Ring Remover with protective cap, "AA" battery and two burrs

K2-4903 chuck for K2-4900

K2-4909 burr, 0.5mm

K2-4910 burr, 1.0mm



Pterygium Burr

K2-4912 Pterygium Burr with protective cap, "AA" battery and 3.5mm burr

K2-4913 replacement burr and chuck assembly for K2-4912



Katena

Ophthalmic Burr
for polishing pterygiums

K2-4920 Includes one 3.3mm diameter round diamond burr, two "AA" batteries and storage case.

K2-4923 diamond burr, round 3.3mm diameter

K2-4925 diamond burr, round 5.0mm diameter

K2-4935 diamond burr, disc shaped 5.0mm diameter

K2-4921 optional lithium hi-energy battery, 3.6 volt, for increased torque



Blade Breakers • Diamond Knives

Castroviejo

Blade Breaker and Holder
with sliding lock

K2-5000



Swiss Model

Blade Breaker and Holder
spherical tip

K2-5150



Diamond LRI Step Knife

6-facet blade for bi-directional cutting
single footplate for better visualization
precalibrated depth settings

K2-6518 Precalibrated
depth setting
of 0.6mm

K2-6519 Precalibrated
depth settings
of 0.5mm, 0.55mm
and 0.6mm

K2-6520 Chu LRI Knife,
0.6mm and
2.5mm preset



K2-6519

600 micron knife also available
sterile disposable K20-2560

Diamond Knives

Diamond LRI Micrometer Knife

6-facet blade for bi-directional cutting
single footplate for better visualization

K2-6522 0 to 1.5mm blade extension
in 10 micron increments



Diamond AK Knife

double
vertical and 45°
cutting edge

K2-6517



Precalibrated depth settings: 0.2, 0.25, 0.3, 0.35, 0.4, 0.6 and 1.5mm



Diamond Micrometer Knife

single
45° cutting edge
1mm wide
200 microns thick

K2-6530



Diamond Micrometer Knife

double
vertical and 45°
cutting edge
200 microns thick

K2-6533



Diamond Micrometer Knife

Thornton "triple"
vertical and 15°
cutting edge with
100 microns base

K2-6535




Diamond Knives

Blade Gauge

For Katena micrometer knives.

K2-6539



K2-6539

Sterilizing Case

with silicone insert
autoclavable plastic


K2-6555 for one knife

K2-6556 for two knives

K9-2024 for three knives



K2-6555



Diamond Knife

30° cutting edge
1mm wide

K2-6560



Diamond Knife

45° cutting edge
1mm wide

K2-6561



Diamond Knife

tri-facet
1mm wide

K2-6565



Diamond Knife

spear
1mm wide

K2-6567



Diamond Knife

10-facet
1mm wide

K2-6568


Diamond Knives



Diamond Step Knife

single 45° cutting edge
1mm wide

K2-6540



Diamond Step Knife

tri-facet
1mm wide

K2-6546



Diamond Step Knife

10-facet
1mm wide

K2-6548

Precalibrated depth settings: 0.2, 0.3, 0.4, 0.5, 0.6, 0.8 and 4.0mm



Diamond Crescent Knife

crescent shaped blade
2.5mm wide, angled

K2-6575



Diamond Phaco Knife

keratome shaped blade

K2-6578 2.8mm straight

K2-6579 2.8mm angled

K2-6580 3.0mm straight

K2-6581 3.0mm angled

K2-6582 3.2mm straight

K2-6583 3.2mm angled



Diamond Knives • LASEK Trephines

Phaco-4 Diamond Knife

keratome shaped blade
4 cutting edges, angled

K2-6603 2.5mm wide

K2-6605 2.8mm wide

K2-6607 3.0mm wide

K2-6609 3.2mm wide

K2-6611 trapezoid
2.7 to 3.2mm



TRAPEZOID



MICS Diamond Knife

keratome shaped blade
4 cutting edges, angled

K2-6622 trapezoid, 1.8 to 2.2mm wide

K2-6654 trapezoid, 2.1 to 2.4mm wide

K2-6680 1.8mm wide

K2-6682 2.2mm wide



TRAPEZOID



Camellin Style

LASEK Trephine

270° semi-sharp edge

90° cutout for flap hinge

K2-7810 8mm diameter

K2-7812 9mm diameter



Camellin Style

LASEK Trephine &

Alcohol Retaining Well

K2-7830 8mm trephine
and 8.5mm well

K2-7832 9mm trephine
and 9.5mm well



Trabeculectomy Punches

Gass

Corneoscleral Punch
1.5mm diameter head
0.75mm deep bite
rotatable tip for convenient positioning

K2-9500



Luntz-Dodick

Trabeculectomy Punch
1.0mm diameter head
with bullet shaped tip
for scleral tunnel insertion
0.5mm deep bite
Cutting head can be preset
in one of 4 positions to suit
the surgeon's preference.

K2-9505 complete


K2-9506 replacement cutting head



Kelly

Descemet's Membrane Punch
1.0mm diameter head
0.75mm deep bite

K2-9510




Trabeculectomy Punches

Fukasaku

Micro Trabeculectomy Punch
0.6mm diameter head
with bullet shaped tip
0.3mm x 0.6mm deep bite

K2-9512



Christensen

Punctum Punch
1mm wide bite
For opening stenotic punctums.

K2-9525

