

Accessories

Marking Tape

autoclavable

For identifying surgical instruments.

K9-1000 red	K9-1050 brown
K9-1010 yellow	K9-1060 blue
K9-1020 orange	K9-1070 green
K9-1030 white	K9-1071 lime green
K9-1035 grey	K9-1080 fluorescent orange
K9-1040 black	K9-1090 fluorescent purple




Vented Tip Guards

autoclavable (pk/100)

For protecting tips of surgical instruments.

K9-1205 assorted colors
K9-1215 white
K9-1225 blue
K9-1235 green
K9-1245 red
K9-1255 yellow
K9-1265 orange



Squeeze Handle

Cleaning Tube

used to flush delicate
mechanized instruments
luer lock connector provided

K9-1350



Sterilizing Cases • Plastic

Sterilizing Case for Diamond Knife

autoclavable plastic
with perforated silicone rack

K2-6555 for one knife

K2-6556 for two knives

Millimeters	Inches
160 x 60 x 25	6 1/4" x 2 1/4" x 1"



Plastic Sterilizing Cases

autoclavable plastic
with silicone finger mat
base, lid and mat are perforated
for steam penetration

Catalog #	Size	# of Instruments	Millimeters	Inches
K9-2018	small, narrow	1-2	195 x 60 x 25	7 5/8" x 2 1/4" x 1"
K9-2020	small	2-3	160 x 70 x 25	6 1/4" x 2 3/4" x 1"
K9-2025	medium	4-6	200 x 110 x 25	7 3/4" x 4 1/4" x 1"
K9-2030	large	10-15	265 x 160 x 25	10 3/8" x 6 3/8" x 1"
K9-2040	large, double level	20-25	265 x 160 x 50	10 3/8" x 6 3/8" x 2"



K9-2020



K9-2040

Sterilizing Cases • Plastic

Sterilizing Case

autoclavable plastic
with perforated silicone rack
base, lid and rack are perforated
for steam penetration

Accommodates 1 instrument.

K9-2022

Millimeters

158 x 70 x 25

Inches

6 ¼" x 2 ¾" x 1"



Sterilizing Case for Knives

autoclavable plastic
with silicone racks
base and lid are perforated
for steam penetration
Accommodates 3 knives.


K9-2024

Millimeters

198 x 70 x 25

Inches

7 ¾" x 2 ¾" x 1"



Sterilizing Cases • Stainless Steel

Katena Sterilizing Case, small
stainless steel
base, lid and mat are perforated
for steam penetration
Accommodates 8 to 12
microsurgical instruments.

K9-2100 complete

K9-2120 silicone mat only


Millimeters	Inches
210 x 125 x 20	8 ¼" x 5" x ¾"



Katena Sterilizing Case, medium
stainless steel
base, lid and mats are perforated
for steam penetration
includes 2 silicone mats, rack for
4 cannulas and a rack for 4 delicate
hooks or knives
Accommodates 20 to 25
microsurgical instruments.

K9-2200


Millimeters	Inches
320 x 220 x 20	12 ½" x 8 ¾" x ¾"



Katena Sterilizing Case, large
stainless steel
base, lid and mats are perforated
for steam penetration
includes 3 silicone mats, racks for 8 cannulas and
a rack for 4 delicate hooks or knives
Accommodates 30 to 40
microsurgical instruments.

K9-2300


Millimeters	Inches
360 x 295 x 20	14 ¼" x 11 ½" x ¾"



Sterilizing Cases • Aluminum

Aluminum Sterilizing Cases, single autoclavable anodized aluminum with silicone finger mat base, lid, and mat are perforated for steam penetration

Catalog #	Size	# of Instruments	Millimeters	Inches
K9-2320	small	6-8	206 x 115 x 19	8 1/8" x 4 1/2" x 3/4"
K9-2330	medium	10-17	247 x 206 x 19	9 3/4" x 8 1/8" x 3/4"
K9-2340	large	16-20	327 x 225 x 25	12 7/8" x 8 7/8" x 1"




Silicone Finger Mat

with "S-Channels" on underside Utilizing a special design which improves air flow beneath the mat, making the drying process more efficient and thorough.

K9-2325 for use with Instrument Case, aluminum, small (K9-2320)

K9-2335 for use with Instrument Case, aluminum, medium (K9-2330)

K9-2345 for use with Instrument Case, aluminum, large (K9-2340)



Sterilizing Cases • Aluminum

Aluminum Sterilizing Case, double anodized aluminum double level with a side compartment for larger items includes removable tray and 4 silicone mats base, lid, tray and mats are perforated for steam penetration Accommodates 35 to 50 microsurgical instruments.
K9-2350

Millimeters	Inches
325 x 265 x 42	12 ¾" x 10 ½" x 1 ½"



Instrument Rack

anodized aluminum Holds and protects delicate hooks, manipulators and cannulas. Fits into upper or lower tray of double level sterilizing case (K9-2350) shown above Accommodates 5 cannulas, 6 hooks, manipulators, etc.
K9-2355

Millimeters	Inches
202 x 115 x 15	8" x 4 ½" x ½"



INSTRUMENT RACK IN TRAY


INSTRUMENT RACK

Eye Shields

Fox

Aluminum
Eye Shield
(bx/50)

K9-6000



¾

Fox


Aluminum Eye Shield
with garter
assorted colors
(bx/50)

K9-6010 adult

K9-6015 pediatric

K9-6030 garter only, assorted*

K9-6040 garter only, white*



¾

*Caution: Garters contain latex which may cause allergic reactions.

Universal

Eye Shield
clear plastic
(pk/12)

K9-6050



¾

Pressure Reducer

Pin Hole
 Eye Shield
 white plastic
 (pk/12)
K9-6100



Honan Balloon
 Intraocular Pressure Reducer
 with relief valve
K9-8410 complete*
K9-8411 headband only, reusable*
K9-8412 headband, latex free
 disposable (bx/12)



*Caution: This product contains latex which may cause allergic reactions.

Tonometer

Schiotz
Tonometer
autoclavable
with certificate of accuracy
K9-8500

