

RETRACTORS





G-16080 *DESMARRES*
LID RETRACTOR
for infants
7 mm wide
overall length 133 mm

02

LID



G-32120 *DESMARRES*
LID RETRACTOR
fig. 0
9 mm wide
overall length 138 mm



G-16060 *DESMARRES*
LID RETRACTOR
fig. 1
12 mm wide
overall length 133 mm



G-16065 *DESMARRES*
LID RETRACTOR
fig. 2
14 mm wide
overall length 135 mm



G-16070 *DESMARRES*
LID RETRACTOR
 fig. 3
 16 mm wide
 overall length 138 mm



G-16075 *DESMARRES*
LID RETRACTOR
 fig. 4
 17 mm wide
 overall length 138 mm



G-32115 *GIESSEN MODEL*
LID RETRACTOR
 14 mm wide, fenestrated
 overall length 133 mm



G-17040 *JAEGER*
LID PLATE
 stainless steel
 20 mm and 23 mm concave ends
 overall length 110 mm



G-15542 *AXENFELD*
LACRIMAL SAC RETRACTOR
4 prongs, blunt
9 mm wide
overall length 140 mm



G-15540 *AXENFELD*
LACRIMAL SAC RETRACTOR
3 prongs, 7 mm wide
overall length 140 mm
blunt
G-15545 sharp



G-15550 *AXENFELD*
LACRIMAL SAC RETRACTOR
2 prongs, 4 mm wide
overall length 140 mm
blunt
G-15555 sharp



G-15560 *ROLLET*
LACRIMAL SAC RETRACTOR
4 prongs, sharp
6 mm wide
overall length 140 mm

02

LACRIMAL SAC



GUTHRIE
RETRACTOR
 double skin hook, sharp
 overall length 125 mm
 small
 large

G-15350
 G-15355



AMERICAN MODEL
LACRIMAL SAC RETRACTOR
 4 x 4 prongs, blunt
 cross-action
 adjustable, self retaining
 overall length 100 mm

G-17085



MUELLER
LACRIMAL SAC RETRACTOR
 3 x 3 prongs
 adjustable
 overall length 70 mm

G-17070



GOLDSTEIN
LACRIMAL SAC RETRACTOR
 3 x 3 prongs
 adjustable
 overall length 30 mm

G-17080